# High-fidelity, high-isotropic resolution diffusion imaging through gSlider acquisition with $B_1^+$ & $T_1$ corrections and multi-coil $B_0$ shim array

Congyu Liao<sup>1,2</sup>, Jason Stockmann<sup>1,2</sup>, Qiyuan Tian<sup>1,2</sup>, Berkin Bilgic<sup>1,2</sup>, Mary Kate Manhard<sup>1,2</sup> Lawrence L Wald<sup>1,2</sup>, and Kawin Setsompop<sup>1,2</sup>

<sup>1</sup> Athinoula A. Martinos Center for Biomedical Imaging, Massachusetts General Hospital, Charlestown, MA, USA.
<sup>2</sup> Department of Radiology, Harvard Medical School, Charlestown, MA, USA

#### Synopsis

gSlider is an SNR-efficient simultaneous multi-slab acquisition that has shown great potential for high-resolution diffusion imaging (DI). In this work, approaches to improve the fidelity of gSlider are proposed. A modified reconstruction which incorporates  $B_1^+$  inhomogeneity and  $T_1$  recovery information was developed and demonstrated to successfully mitigate slab-boundary artifacts in short-TR gSlider acquisitions. Slice-by-slice  $B_0$ -shimming through multi-coil shim-array, and high in-plane acceleration through virtual-coil GRAPPA were also incorporated into the acquisition and demonstrated to achieve 8-11x reduction in  $B_0$  distortion in single-shot EPI. The modified gSlider acquisition/reconstruction was used to acquire high-fidelity whole-brain 1mm isotropic DI with 64 diffusion-directions in 20-minutes at 3T.

#### Introduction

High resolution diffusion imaging (DI) is limited by low-SNR. Generalized slice dithered enhanced resolution (gSlider) is an SNR-efficient acquisition technique with self-navigated RF slab-encoding, which has shown great potential for high resolution DI (1). In this work, we developed approaches to improve the fidelity of gSlider; to make it robust to slab-boundary artifacts, and to drastically reduce distortion from B<sub>0</sub> inhomogeneity. To mitigate slab-boundary artifacts in short-TR acquisitions, a dictionary-based B<sub>1</sub><sup>+</sup>-inhomogeneity and T<sub>1</sub>-recovery correction is incorporated into gSlider reconstruction. To mitigate image distortions, a 32-channel integrated  $\Delta$ B0/Rx array (AC/DC coil) (2) is utilized for slice-by-slice shimming, to reduce B<sub>0</sub>-inhomogeneity by >50%. This is then combined with high in-plane acceleration (R<sub>inplane</sub>=4) and virtual conjugate coil (VCC-) GRAPPA reconstruction(3,4), to achieve 8-11x distortion reduction in single-shot gSlider-EPI. We demonstrated that the proposed method can achieve high-fidelity whole-brain 1mm-isotropic DI with 64 diffusion-directions in **20 minutes** on a 3T scanner.

#### Methods

Figure1(a) shows the sequence diagram of gSlider, where two external triggers are added in each TR to enable slice-by-slice  $B_0$ -shimming with the AC/DC coil (Fig1(b)). To avoid poor performance in whole-brain fat suppression from large out-of-slice  $B_0$ -inhomogeneity, the slice-by-slice shimming is turned off during fat saturation. In each TR, an additional  $k_y$  blip is

added to shift k-space and create more unique source points for improved VCC-GRAPPA reconstruction (4). To implement slice-optimized shimming, a field-map with conventional global-shim is acquired. Based on shim optimization on this field-map, the shim-arrays generate slice-specific DC currents for dynamic updating. For simultaneous multi-slice acquisition, the simultaneously acquired slices are shimmed jointly.

To mitigate slab-boundary artifacts, RF-encoding imperfections due to  $B_1^+$ -inhomogeneity and incomplete  $T_1$ -recovery are estimated and incorporated into the RF-encoding matrix (forward model) of the gSlider reconstruction. Figure2(a) shows the flowchart of  $B_1^+$ -inhomogeneity correction, where RF-encoding profiles at a range of flip-angles are Bloch simulated and assigned to different spatial locations represented within the RF-encoding matrix (5) using a discretized  $B_1^+$  map. For incomplete  $T_1$ -recovery, non-ideal slab-profiles of RF-encodings can cause partial excitations in adjacent slabs which are not fully-recovered in slab-interleaved acquisitions with short-TRs. This effect is also modeled into the Bloch simulation of the RF-encodings (assuming average  $T_1=1000$ ms), and incorporated into the encoding matrix. Figure 2(c) shows the partial  $M_z$ -recovery from adjacent slabs excitations before and after TR/2 longitude-relaxation at various  $B_1^+$  excitation levels.

<u>Acquisition</u>: the following data were acquired using a Siemens Prisma 3T-Scanner with a custom AC/DC array:

i) To validate slice-optimized shimming, gSlider-EPI with 5 slab-encodings were acquired: FOV:  $220 \times 220 \times 160 \text{ mm}^3$ , 32 thin-slabs (5mm slab-encoding), TR/TE=5100/77ms, echo spacing=0.93ms. To assess distortion, data were acquired using both Anterior-to-Posterior (AP) and Posterior-to-Anterior (PA) phase-encodings at different R<sub>inplane</sub> accelerations, with and without slice-optimized shimming. A matching T<sub>2</sub>-SPACE data was also acquired as a distortion-free reference.

ii) whole-brain gSlider-EPI with dynamic shimming: 1mm isotropic with 64 diffusion-directions,  $MB \times R_{inplane} \times gSlider = 2x4x5$ , FOV:220x220x160mm<sup>3</sup>, b=1000s/mm<sup>2</sup>, TR/TE=3500/77ms,  $T_{acq}=20min$ .

 $B_1^+$  and  $B_0$  maps were obtained at 2mm in-plane resolution.  $B_1^+$  map was acquired using Turbo-FLASH with 4mm slice-thickness and 25% gap.  $B_0$  map was acquired using multi-echo GRE with 2.5mm slice-thickness and 100% gap. The slice resolution including gap for both matches the 5mm gSlider slab-encoding.

<u>Reconstruction and post-processing</u>: VCC-GRAPPA with phase-matching (4) was used for parallel imaging reconstruction. The 5 RF-encoded volumes of each diffusion-direction were then combined to create 1mm isotropic data, using gSlider reconstruction with proposed modified RF-encoding matrix. The data were then motion&eddy-current corrected, and averaged-DWI and FA-maps generated using FSL toolbox(6).

## Results

Figure2(c) shows the  $B_1^+$  and  $T_1$  corrections in gSlider reconstruction, where slab-boundary artifacts are well mitigated for a short-TR acquisition at 3.5s. Figure3 show comparisons of the image distortion, with and without dynamic-shimming. The green arrows highlight the distortions that were alleviated, with slice-by-slice shimming achieving >50% reduction in  $\Delta B_0$ . This shimming is then combined with  $R_{inplane}4$  to achieve 8-11x reduction in  $\Delta B_0$  distortion, to create images with outlines closely matching that of  $T_2$ -SPACE. Figure4(a&b) shows the  $B_0$ -field from MB=2 dynamic shimming, where two *distant* slices are shimmed jointly, providing similar performance to the MB=1 case. Figure 5 shows averaged-DWIs and FA-maps of the 1mm isotropic DI data. High-quality results are shown, with minimal distortions even in the typically problematic areas of temporal and frontal lobes (Figure 5(c)).

## Discussion and conclusion

In this work, we developed approaches to improve gSlider acquisition where i) reconstruction with  $B_1^+$  and  $T_1$  corrections effectively mitigate slab-boundary artifacts in short-TR acquisitions, and ii) slice-by-slice shimming and high in-plane acceleration achieve 8-11x reduction in  $B_0$ -distortion. The proposed corrections and local-field control were used to acquire high-quality, high-fidelity DI at 1mm isotropic in 20 minutes.

## Acknowledgement

This work was supported in part by NIH research grants: R01 MH116173, R01EB020613, U01EB025162, P41EB015896, and the shared instrumentation grants: S10RR023401, S10RR019307, S10RR019254, S10RR023043

## References

1. Setsompop K, Fan Q, Stockmann J, et al. High-resolution in vivo diffusion imaging of the human brain with generalized slice dithered enhanced resolution: Simultaneous multislice (gSlider-SMS). Magn. Reson. Med. 2018;79:141–151. doi: 10.1002/mrm.26653.

2. Stockmann JP, Witzel T, Keil B, Polimeni JR, Mareyam A, Lapierre C, Setsompop K, Wald LL. A 32-channel combined RF and BOshim array for 3T brain imaging. Magn. Reson. Med. 2016;75:441–451. doi: 10.1002/mrm.25587.

3. Blaimer M, Gutberlet M, Kellman P, Breuer FA, Köstler H, Griswold MA. Virtual coil concept for improved parallel MRI employing conjugate symmetric signals. Magn. Reson. Med. 2009;61:93–102. doi: 10.1002/mrm.21652.

4. Liao C, Manhard MK, Bilgic B, et al. Joint Virtual Coil Reconstruction with Background Phase Matching for Highly Accelerated Diffusion Echo-Planar Imaging. In: ISMRM . Paris; 2018. p. 0465.

5. Pauly J, Le Roux P, Nishimura D, Macovski A. Parameter relations for the Shinnar-Le Roux selective excitation pulse design algorithm (NMR imaging). IEEE Trans. Med. Imaging 1991;10:53–65. doi: 10.1109/42.75611.

6. Jenkinson M, Beckmann CF, Behrens TEJ, Woolrich MW, Smith SM. Fsl. Neuroimage 2012;62:782–790. doi: 10.1016/j.neuroimage.2011.09.015.



Figure 1. (a) sequence diagram of gSlider acquisition and (b) the 32-channel combined RF and BO shim array.



Figure 2. The flowchart of (a)  $B1^+$  inhomogeneity and (b) T1 corrections for gSlider data. (c) the results with and without correction.



Figure 3. The comparisons of the field maps and EPI distortions with Anterior-to-Posterior (AP) and Posterior-to-Anterior (PA) phase-encodings, with and without dynamic shimming. The shimming is then combined with  $R_{inplane}4$  acceleration to achieve ~10-fold reduction in EPI distortion, and create low-distortion images. The red outlines shows the contour of T2-SPACE reference.



Figure 4. The shim results for a SMS slice group. (a) compared to global shim, in MB=2 case, the two distant slices are shimmed jointly, while out-of-slices regions are unconstrained and are allowed to have poor shim. The MB=2 shimming performance closely resembles that of the MB=1 slice-optimized shimming, as shown in (b).



Figure 5. (a) the averaged-DWIs and (b, c) FA maps of the 1 mm isotropic DI data with 64 directions obtained in just **~20-minutes** using the proposed correction and dynamic shimming framework. (c) The typically problematic areas of the temporal and frontal lobes of the brain can be observed by using dynamic shimming.