Development of Machine Learning Classifiers for Blood-based Diagnosis and Prognosis of Suspected Acute Infections and Sepsis

Ljubomir Buturović* **Michael Mayhew Roland Luethy** Kirindi Choi Uroš Midić Nandita Damaraju Yehudit Hasin-Brumshtein **Amitesh Pratap** Rhys M. Adams Joao Fonseca **Ambika Srinath** Paul Fleming **Claudia Pereira Oliver Liesenfeld Purvesh Khatri Timothy Sweeney**

Abstract

We applied machine learning to the unmet medical need of rapid and accurate diagnosis and prognosis of acute infections and sepsis in emergency departments. Our solution consists of a MyrnaTM Instrument and embedded TriVerityTM classifiers. The instrument measures abundances of 29 messenger RNAs in patient's blood, subsequently used as features for machine learning. The classifiers convert the input features to an intuitive test report comprising the separate likelihoods of (1) a bacterial infection (2) a viral infection, and (3) severity (need for Intensive Care Unit-level care). In internal validation, the system achieved AUROC = 0.83 on the three-class disease diagnosis (bacterial, viral, or non-infected) and AUROC = 0.77 on binary prognosis of disease severity. The Myrna, TriVerity system was granted breakthrough device designation by the United States Food and Drug Administration (FDA). This engineering manuscript teaches the standard and novel machine learning methods used to translate an academic research concept to a clinical product aimed at improving patient care, and discusses lessons learned.

1 Introduction

New advances on research in applications of machine learning (ML) and artificial intelligence in medicine are published on a regular basis. However, there is lack of literature on translation of these innovations to clinical practice, in particular for tests based on molecular data. To fill the gap, we report on the development of classifiers for detecting type and severity of infections in patients who present in emergency departments (ED) with symptoms of acute infection and sepsis, an unmet medical need [1].

Current modalities for diagnosing infections mostly rely on detection and identification of pathogens. However, this is inadequate because in the majority of cases, pathogens are not found in blood or anywhere else in the body [2]. A more recent approach relies on response of the immune system to the infection (*host response*). This could provide diagnostic and prognostic information regardless of whether a pathogen is eventually identified.

From the ML perspective, diagnosing infections in the ED can be reduced to two classification problems: infection type classification (diagnosis) and illness severity classification (prognosis). The

^{*}lbuturovic@inflammatix.com. Authors are at Inflammatix Inc., Sunnyvale, CA, United States.

classifiers use gene expression (abundance of mRNAs) of cells from whole blood as input features. The diagnostic classifier estimates the probability of the patient having bacterial, viral or no infection (BVN), whereas the prognostic binary classifier estimates probability of severe outcome in the given time window (SEV).

We used clinical adjudication [3] as the ground truth for the BVN classifier, and 30-day survival as the ground truth for the SEV classifier. We used 29 mRNAs as input numerical features, measured using a variety of measurement platforms during training and validation, comprising both commercially-available (microarrays, RNA-Seq and molecular barcoding technology (NanoString[®]), which are established technologies for measuring gene expression) and our in-development platform Myrna[™], which uses a rapid method called Loop-Mediated Isothermal Amplification (LAMP). The LAMP experiments were performed using two approaches: 1. "benchtop LAMP", which is not fully-automated and uses a commercially available instrument. It was utilized before TriVerity[™] became available, and 2. "automated LAMP" using TriVerity.

The final test produces three scores (Bacterial, Viral, Severity) based on underlying classifier probabilities. The probabilities are divided into five bands, corresponding to very low, low, moderate, high, and very high probabilities of the corresponding class. The probabilities and bands are further transformed to scores on 0-50 scale and presented to the user (Fig 1). The test report was developed using in-house expertise, feedback from ED and clinical laboratory physicians, FDA predicate device clearances, and additional market research.



Figure 1: A sample TriVerity test report. The scores are monotonic piecewise-linear functions of the probabilistic outputs of the BVN and SEV classifiers. During the development, severe illness was defined as death within 30 days of the test administration.

The system software overview is described in Section 8.1.

In this paper, we describe how the classifiers were developed and their performance characteristics. Inspired by [4], we show the key elements of the classifier development in Fig. 2. As in [4], we note that the classifier selection and tuning represent a modest fraction of the overall effort from academic prototype to product launch. The manuscript sections correspond to the elements shown, except where noted. Given the extensive validation of the 29-mRNA signature in academic prototypes [5, 6, 7], here we specifically emphasize activities specific to commercially usable point-of-care classifier development.

2 Development data and classifier development

First we describe the development data, with emphasis on properties that have significant impact on the ML methods: batch effects and platform transfer.

2.1 Data gathering and adjudication

The sources of development data were public (generated using microarrays, from GEO [8] and ArrayExpress [9]) and proprietary (generated using NanoString[®] and qPCR) [10]. The approach



Figure 2: **Key elements of the clinical classifier development process.** The white rectangles are generic machine learning activities. The blue rectangles are specific to clinical classifiers. The green rectangles are specific to clinical classifiers based on molecular features. QSR: Quality Systems Regulation; FDA: United States Food and Drug Administration. The names correspond to the manuscript sections, except *Clinical endpoints* which will be described elsewhere.

to selecting the development data was guided by prior research on diagnosing infections using a molecular approach [6].

The development data was split into training and validation sets. The training set was used for hyperparameter tuning with cross-validation (CV), while the validation data was used for comparison of the tuned classifiers, and never for classifier training nor hyperparameter tuning (except exactly once to train the final locked classifiers). This approach reduces bias compared with CV [10].

The class labels (*clinical truth*) for the development data were derived from the public records for public data and from a clinical adjudication process for proprietary data. As such, the public data were inherently more noisy than the proprietary data, but also better representative of the heterogeneity in the real-world patient population. The proprietary data were from company-sponsored clinical studies, all with physician-adjudicated label assignment. The proprietary studies were approved by institutional review boards.

The detailed information about the development datasets is given in Section 8.2.

2.2 Batch effects

The blood samples used to generate data for the development were processed in different hospitals across independent studies and generated using different technical platforms, on biologically and clinically heterogenous patient populations (ED, intensive care units, different races and ethnicities, different disease severity). We used normalization to adjust different levels of gene expression, but the normalization does not correct for technology- and study-induced batch effects. To reduce their impact, we used validation and CV based on *grouped* approach.

To reduce the platform-induced batch effects while maintaining biological and clinical heterogeneity, we grouped the training and validation data by gene expression platform. The training data were almost exclusively ($\approx 98\%$) microarray data, whereas the validation data were assayed on NanoString (during development) or the target platform (for production classifiers). The goal was to ensure the robustness of the trained classifiers across a variety of platforms. The advantage of this robustness was that it allowed use of large amounts of publicly available datasets as they become available to further improve accuracy of our classifiers, while maintaining our ability to translate them to a point-of-care platform. The gap between the training and validation performance served as an indicator of the classifier's robustness with respect to the gene expression platform.

To reduce the study-induced batch effects, we used grouped CV [11]. The grouped CV was used in the hyperparameter tuning to evaluate performance of candidate hyperparameter configurations (HC). In grouped CV, the folds were configured such that different folds contained disjoint sets of studies (in other words, all blood samples from a given study were assigned to the same grouped CV fold). Analogous to the approach used for platform-induced batch effects, HCs which exhibited poor performance on the left-out folds (usually including study batch effects) were automatically eliminated by the hyperparameter tuning algorithms.



Figure 3: **Overview of classifier tuning.** CF: Concordance Filtering (see Section 2.3); HC: Hyperparameter Configuration.

2.3 Hyperparameter tuning, classifier selection and platform transfer

An overview of hyperparameter tuning is shown in Fig. 3.

A key requirement of our classifiers was *platform transfer*, defined as generalization from the mostly public expression data, which comprised the training set, to the target instrument. During most of the classifier development, the target instrument was not available because the device itself was being developed: the product development took about six years from conception to the design freeze (a timepoint at which device is considered finalized), whereas the instrument data first became available about six months prior to the design freeze. As a result, input features for the classifier training were measured using different, commercially available instruments (NanoString[®] was used as a proxy for MyrnaTM during the development). The platform transfer problem can be defined as ensuring that the classifiers developed on the available platforms would perform similarly on the target platform, once it becomes available. To enable the platform transfer, we developed a *Concordance Filtering* (CF) approach.

The platform transfer based on CF supports an efficient use of the heterogenous data. It is an enabling step for applying the classifiers developed on one or more platforms to a different or newly developed platform. The number of samples available on our target platform was relatively limited and deemed insufficient alone for training the classifiers that would generalize to a global population. The goals of the platform transfer were A) use all available data, on all platforms, for classifier training, and B) maximize performance and generalizability of the classifier on the target platform.

To develop an approach which satisfies the goals, we relied on the following ideas:

- In early versions of the classifiers, optimize for adequate performance on the NanoString data as the highest-quality feature values available
- Later, jointly optimize for high level of concordance (i.e., similar gene expression values measured across replicate samples) between NanoString and Myrna
- Thus, tuning the classifiers to achieve good concordance between platforms can be performed using relatively small unlabeled sets of samples assayed on both platforms

The CF algorithm removes, during hyperparameter tuning, the classifiers which did not achieve adequate concordance between the replicate samples. The concordance was measured by the coefficient of determination (R^2) and/or Pearson correlation coefficient [PCC] between NanoString[®]- and MyrnaTM-generated probabilities. In other words, during the tuning, we applied a classifier to the set

of replicate NanoString/Myrna samples. If the concordance was below a fixed minimal value, the classifier was removed from the tuning process, regardless of clinical performance. For BVN, where there are three classes, the CF filtering step was based on the minimum PCC or R^2 for bacterial or viral class, which are of primary interest because they are reported to the user. For SEV, the CF filtering was based on the probability of the severe class.

Using this approach, we developed the instrument and classifiers in parallel, thereby significantly reducing the product development timelines and cost, and substantially boosting the number of samples available for training. The final tuning of the classifiers using Myrna data was done in the very last stages of the project and required a modest number of samples. Specifically, for the last versions of the classifiers we assayed a set of matching (replicate) blood samples on NanoString and Myrna[™] for use in the CF.

2.4 Reproducibility

Reproducibility refers to an important safety property of diagnostic tests to generate clinically similar results using repeated measurements from the same patient (for example, multiple tubes from a single blood draw). It is conceptually similar to platform transfer except that the replicates come from repeated measurements on the same platform. Accordingly, we adopted the following approach: during CF, add Myrna replicates to the set of NanoString/Myrna replicates. That way, we ensured that the highly-ranked classifiers achieved high level of reproducibility within and across platforms.

2.5 Hyperparameter tuning: ensembling and loss

Our approach used ensembling. We evaluated several variants of ensembling to empirically find the best performing method. The hyperparameter tuning was repeated M = 10 times independently, with different random seeds and otherwise identical configurations. The class probabilities were the averages of the M output probabilities generated by the ensemble members.

The loss function used for ranking the HCs was Area Under Receiver Operating Characteristic (AUROC) for severity, and the multi-class generalization of AUROC [12] (mAUROC) for BVN classifier. Ideally, the candidate classifiers should be ranked using the performance statistics of the final product, such as sensitivity and specificity. However, the clinical statistics (Section 3) involve setting of decision thresholds, which is not amenable to full automation.

During CV, the mAUROC can be computed using averaging or pooling. In averaging, mAUROC is estimated for each fold, and the final mAUROC is the average across folds. In pooling, predicted probabilities for each fold are recorded, and the final mAUROC is computed by pooling the fold probabilities and estimating mAUROC using the pooled data. While averaging may be preferred [13, 14], in our case it was not an option because the majority of studies did not have patients in all three classes.

Based on expert judgment, we considered an increase in classification performance of 0.01 AUROC units to be clinically significant.

We used Optuna [15] to select candidate hyperparameter configurations (HC) from a search space of 1000. We chose a fixed value of five folds for all CV experiments.

2.6 Assay-guided feature selection

Academic research identified 29 candidate genes (a *signature*) for diagnosis and prognosis of presence, type and severity of infectious diseases [5, 6, 7]. They formed an initial set for the search for genes compatible with the target device. The search is described in Section 8.3.

2.7 Classifier lock

The classifier lock is the setting of classifiers' parameters (weights) to final optimized values for deployment in clinical care. A key practical decision was to choose between the following options:

1. Pool the training and validation data and train the final classifier using the pooled dataset with hyperparameters found by hyperparameter tuning. In this case, training data is maximized, and includes the target platform data, but there is no data left to validate the locked classifiers.

2. Lock the models trained on training data. The classifiers are validated using the validation set. However, the training data is reduced, and does not use the target platform.

We opted for pooling, based on the strong preference to incorporate the target platform in the training of the locked classifiers.

2.8 Fairness analysis

The performance of the development versions of the classifiers across race, age and sex subgroups was reported in [16]. The performance of the production classifiers will be published in a future manuscript, following the completion of the pivotal trials.

2.9 QSR/FDA compliance

The Myrna[™] system and TriVerity[™] Test, including the classifiers, are subject to the Quality Systems Regulation by the United States FDA. The regulation is codified in Title 21 of the Code of Federal Regulations (CFR) Part 820 [17] and is meant to ensure quality and safety of medical devices intended for human use. It is a complex topic outside of scope of this manuscript. Machine learning activities impacted by the Part 820 include documentation, personnel training, corrective and preventive actions, hazard analysis, statistical techniques and classifier updates.

3 Clinical performance statistics

The key test performance characteristics comprise likelihood ratios for individual test bands, band sensitivities and specificities, coverage and monotonicity.

For tests with two possible test results, positive and negative (binary tests), the positive diagnostic likelihood ratio (LR) is defined as the ratio of the probability of a positive test result in patients with the disease to the probability of a positive test result in patients without the disease. The negative diagnostic ratio is defined as the ratio of the probability of a negative test result in patients without the disease. The negative disease to the probability of a negative test result in patients without the disease [18]. The two LRs are the properties of a classifier and a decision threshold, and fully characterize the test performance. However, TriVerity has five possible test results (assigned bands, Fig 1). This impacts the definitions and interpretations of the standard performance statistics. The definition of the LRs in such cases is non-obvious because it is not apparent what constitutes positive and negative test results for the bands. For computing the LRs for band i, we considered all samples assigned by the classifiers to band i as *predicted positive*, and samples assigned outside the band as *predicted negative*. This permits an intuitive behavior of monotonically increasing likelihood ratios from band 1 to band 5 (an established desirable property of multi-band tests [19]). Given the above convention, the positive LRs become:

$$LR_{i} = \frac{\frac{x_{i}}{X}}{\frac{y_{i}}{Y}}, i \in 1..5, X = \sum_{i=1}^{5} x_{i}, Y = \sum_{i=1}^{5} y_{i}$$
(1)

where x_i, y_i are the number of samples with ground truth **POS** and **NEG** in band *i*, respectively.

4 Decision threshold tuning

Decision thresholds help clinicians make diagnosis and treatment decisions. Each TriVerity output probability is split in five bands using four thresholds (Fig 1).

The decision threshold tuning required the specification of the threshold-tuning-set (TTS) and a decision tuning algorithm. The assumption is that the samples in TTS have been analyzed using a classifier, producing predictive probabilities. Subsequently, we tune decision thresholds using the TTS predicted probabilities and the combination of decision tuning algorithm and manual finetuning.

In problems with two classes and a binary test report, the threshold tuning is typically achieved by finding a clinically best trade-off between sensitivity and specificity, for example by expert judgment. However, TriVerityTM requires setting of 8 thresholds for BVN and 4 for SEV (Fig 1). To facilitate the

threshold selection during model tuning, we introduced an algorithm called GAT (Genetic Algorithm for Thresholds, [20]). GAT is based on genetic (evolutionary) optimization as implemented in the DEAP python library [21]. Each individual (*chromosome*) in GAT is a set of the decision threshold values. GAT operates by maximizing a fitness function through evolutionary computation until the chromosome with the best fitness is found. GAT is applied to each test output (bacterial, viral, severity) independently. It is defined as follows:

Input: desired (target) values for LR1, LR5, coverage, and percent moderate

Output: best thresholds (4 values)

- 1. The initial population for the evolutionary (genetic) algorithm is randomly generated. A set of chromosomes, each representing a potential solution to the problem, is created. The chromosome corresponding to a solution has four values, representing the four thresholds needed to split the probabilities into five bands.
- 2. The fitness of each chromosome in the population is evaluated using a fitness function. The function assigns a fitness score to each chromosome based on how well it fits the desired criteria for LR1, LR5 and coverage. The coverage is only considered if the value is below the target value, so that coverage exceeding the target is not penalized.
- 3. A new generation is created by selecting parents according to their fitness. Offspring are created using crossover and mutation operations. The individuals with the top 20% fitness are always kept in the population.
- 4. Steps 2 and 3 are repeated for a given number of iterations.

The specification of TTS is, to the best of our knowledge, an unsolved problem. The problem is that, following classifier development, a separate dataset may be required to tune the decision thresholds, and that may not always be available. We considered the following options: A) using the same training data for threshold tuning and performance evaluation and B) setting aside a separate set of data exclusively for TTS. The approach A is a resubstitution approach, which may introduce optimistic bias. The approach B may not use the valuable data efficiently. We decided between the two approaches as follows:

- Tune the thresholds using resubstitution and GAT algorithm. This produces thresholds T
- Estimate clinical performance R using TTS probabilities and thresholds T. R is a set of relevant statistics such as LRs and coverage
- Tune the thresholds using a reduced-bias approach (the reduction refers to bias in comparison with the resubstitution approach). This means using different set of data for threshold tuning vs. classifier performance assessment
- Record the resulting clinical performance U
- If the difference between R and U is clinically insignificant, use the data-efficient resubstitution approach for decision thresholds tuning.

The algorithm requires reduced-bias threshold selection. To that end, we introduced the following algorithm for decision threshold CV:

- Split TTS (i.e., probabilities) into K folds. Each fold consists of learning set (LSET) and cross-validation set (CVSET)
- For each (LSET, CVSET) pair:
 - Tune decision thresholds using GAT algorithm and LSET dataset
 - Assign the CVSET probabilities to different bands using the tuned thresholds
- Compute the TTS clinical performance using the pooled band assignments

The threshold CV computes reduced-bias estimate of the clinical performance when thresholds are tuned using independent data. This estimate is then used to decide whether the data-efficient resubstitution method has acceptable bias. If the conclusion is positive, we can use the same dataset for threshold tuning and estimation of performance; otherwise, a separate TTS is needed.

Table 1: **Classifier accuracy and platform-to-platform concordance comparison.** MLP: Multilayer Perceptron; LOGR: Logistic Regression; RBF: Support Vector Machine with Radial Basis Function kernel; SVM: Support Vector Machine with linear kernel; LGBM: Light GBM classifier; XGB: XGBoost; GPC: Gaussian Process Classifier. Concordance: Pearson correlation coefficient between replicate samples assayed on NanoString[®] and LAMP. The classifiers are sorted by validation mAUROC. The ties are resolved by concordance.

Algorithm	CV mAUROC [95% CI]	Val. mAUROC [95% CI]	Concordance [95% CI]
MLP	0.88 [0.88, 0.89]	0.84 [0.81, 0.86]	0.94 [0.93, 0.94]
LOGR	0.89 [0.88, 0.89]	0.84 [0.81, 0.86]	0.93 [0.92, 0.94]
RBF	0.87 [0.86, 0.88]	0.83 [0.80, 0.86]	0.95 [0.94, 0.96]
SVM	0.87 [0.86, 0.87]	0.83 [0.80, 0.86]	0.92 [0.91, 0.93]
LGBM	0.87 [0.86, 0.88]	0.83 [0.80, 0.86]	0.88 [0.86, 0.89]
XGB	0.87 [0.86, 0.88]	0.82 [0.78, 0.84]	0.80 [0.78, 0.82]
GPC	0.86 [0.85, 0.87]	0.81 [0.77, 0.84]	0.76 [0.73, 0.78]

The selection of thresholds during classifier development used GAT due to good performance, automation and high throughput; however, for the final classifier lock, the thresholds were manually finetuned, to achieve the best trade-off between different clinical performance statistics.

5 Results

The quality of input data and normalization approaches were critical to the development. The BVN mAUROC improved from under 0.7 in the early days of the development, which was not viable for launching the product, to 0.83, which provides favorable overall clinical performance. We hypothesize that the progress was achieved by increasing the quality and quantity of the public and proprietary data, and improving development methodology. The data increase qualitatively changed the type of top-ranked classifiers, which oscillated between multi-layer perceptron (MLP) and LOGR, until the introduction of platform transfer. At design freeze, the BVN development data comprised 4306 training and 679 validation samples, derived from 47 and 7 studies, respectively. The SEV development data comprised 2758 training and 723 validation samples, derived from 38 and 8 studies, respectively.

We periodically evaluated seven leading classification algorithms at different time points during the development. A late-stage snapshot is shown in Table 1. The winning classifier was LOGR for both BVN and SEV problems.

The greedy forward search did not yield a subset of genes which improved clinical or analytical performance for at least one of the classifiers (results not shown). Based on this observation, we decided to finalize the product using 29 genes.

The final results of TriVerityTM classifiers are shown in Table 2. Importantly, we found that for N = 534 patients for which both BVN and SEV ground truths were available, 73% were assigned by the classifiers to at least one of the outer bands (1 or 5). This is clinically important because those diagnoses provide the most confidence in their result. In contrast, only 0.2% of the patients in the common set were assigned to band 3 for all three scores. This is important because the band 3 has LRs that are clinically uniformative.

We observed 0.886 BVN mAUROC in training and 0.833 in validation. We hypothesize that the principal causes of the gap are the major platform difference (BVN training is almost entirely microarray data, whereas validation is MyrnaTM data) and significant differences in patient populations.

We observed 0.886 SEV AUROC in training and 0.793 in validation. Besides platform gap, the severity gap was due to inclusion of pediatric patients in severity training set. The rationale for this design choice is that removal of pediatric patients from training had limited to no impact on the validation results, and it may improve the performance in a future expanded indication for the test.

To estimate the impact of CF, we performed an ablation study in Section 8.4.

Table 2: **Clinical performance of the winning TriVerity**TM **classifiers.** All samples were processed using MyrnaTM Instrument. PPA: Positive Percent Agreement (equivalent to sensitivity when ground truth is not 100% reliable); NPA: Negative Percent Agreement (equivalent to specificity when ground truth is not 100% reliable). Coverage is the percent of patients assigned to the indicated band(s). *Moderate* refers to band 3.

Band(s)	Statistic	Training Value [95% CI]	Val. Value [95% CI]
	BVN mAUROC	0.886 [0.877, 0.894]	0.833 [0.806, 0.861]
1	Bacterial LR	0.09 [0.07, 0.11]	0.08 [0.03, 0.14]
5	Bacterial LR	7.1 [6.2, 8.2]	7.6 [4.8, 15]
1	Bacterial PPA	95 [94, 97]	98 [96, 99]
5	Bacterial NPA	93 [92, 94]	95 [93, 98]
1, 5	Bacterial coverage	60 [58, 61]	35 [31, 39]
3	Bacterial moderate	7.1 [6.3, 7.8]	16 [13, 19]
1	Viral LR	0.07 [0.05, 0.10]	0.13 [0.05, 0.22]
5	Viral LR	8.1 [7.3, 9.1]	9.6 [7.0, 14]
1	Viral PPA	97 [96, 98]	95 [92, 98]
5	Viral NPA	90 [89, 91]	93 [91, 95]
1, 5	Viral coverage	54 [53, 56]	52 [48, 55]
3	Viral moderate	13 [12, 14]	15 [12, 18]
1, 5	BV coverage	72 [71, 74]	58 [54, 62]
	SEV AUROC	0.886 [0.866, 0.906]	0.793 [0.738, 0.843]
1	SEV LR	0.04 [0.01, 0.09]	0.16 [0.04, 0.33]
5	SEV LR	7.3 [5.8, 9.0]	8.9 [5.2, 16]
1	SEV PPA	98 [95, 100]	94 [89, 99]
5	SEV NPA	93 [92, 94]	97 [95, 98]
5	SEV coverage	61 [59, 63]	40 [36, 44]
3	SEV moderate	6.5 [5.5, 7.4]	15 [12, 17]

6 Discussion

We performed successful transfer of a research prototype of a clinical infectious disease diagnostic test to a design freeze product ready for independent verification and validation, across heterogenous data environment. The results suggest that the BVN and SEV ML classifiers perform adequately for clinical deployment [22]. Approximately three out of four patients are likely to receive at least one clinically highly informative score (bands 1 or 5), and almost none (0.2% in the validation set) are likely to receive an uninformative test report. Based on our assessment of the clinical utility and commercial potential of the validation results, we built and locked the production classifiers and initiated pivotal trials.

We found that concordance filtering significantly improved the quality of platform transfer for some classification algorithms (MLP). For others, like LOGR, it served as assurance that the winning classifier would meet the platform transfer requirements.

The validation performance showed that grouped CV effectively countered the platform and study batch effects. While the effects were not completely eliminated, they were nevertheless sufficiently attenuated to enable clinical application of the system.

The classifier comparison showed relatively small but clinically relevant differences in performance by AUROC. The gradient-boosting classifier (XGBoost and Light GBM), generally considered top choice for classification of tabular data, underperformed in our use case, both by classification performance and concordance among platforms. This is a somewhat surprising finding, potentially caused by strictly numeric type of the input features. We also observed consistent underperformance by Gaussian Process Classifier, which may be useful learning because this classifier is relatively rarely evaluated in ML manuscripts. These findings are under active investigations.

Arguably, our use of pooling the predicted probabilities across CV folds is biased. However, the clinical application and the composition of development data necessitated this approach. We found the resulting performance to be robust in validation set and across gene expression platforms.

Our work and manuscript have limitations. We presented development data; a prospective evaluation of the clinical validity is currently in progress. The clinical utility of TriVerity[™] will be determined in prospective and randomized trials, currently in planning stages. The adjudicated ground truth for TriVerity development is imperfect and may affect the performance statistics reported here.

7 Conclusion

We presented the development of clinical classifiers for rapid diagnosis and prognosis of acute infections and sepsis in emergency departments based on blood draw. Internal evaluations suggest that the product is adequate for clinical use and commercial launch, pending clearance by the United States FDA.

References

- [1] W. Bauer, K. Kappert, N. Galtung, D. Lehmann, J. Wacker, H. K. Cheng, O. Liesenfeld, L. Buturovic, R. Luethy, T. E. Sweeney *et al.*, "A novel 29-messenger rna host-response assay from whole blood accurately identifies bacterial and viral infections in patients presenting to the emergency department with suspected infections: a prospective observational study," *Critical care medicine*, vol. 49, no. 10, pp. 1664–1673, 2021.
- [2] S. Jain, W. H. Self, R. G. Wunderink, S. Fakhran, R. Balk, A. M. Bramley, C. Reed, C. G. Grijalva, E. J. Anderson, D. M. Courtney *et al.*, "Community-acquired pneumonia requiring hospitalization among us adults," *New England Journal of Medicine*, vol. 373, no. 5, pp. 415–427, 2015.
- [3] N. N. Whitfield, C. A. Hogan, J. Chenoweth, J. Hansen, E. B. Hsu, R. Humphries, E. Mann, L. May, E. A. Michelson, R. Rothman *et al.*, "A standardized protocol using clinical adjudication to define true infection status in patients presenting to the emergency department with suspected infections and/or sepsis." *Diagnostic Microbiology and Infectious Disease*, p. 116382, 2024.
- [4] D. Sculley, G. Holt, D. Golovin, E. Davydov, T. Phillips, D. Ebner, V. Chaudhary, M. Young, J.-F. Crespo, and D. Dennison, "Hidden technical debt in machine learning systems," Advances in neural information processing systems, vol. 28, 2015.
- [5] T. E. Sweeney, A. Shidham, H. R. Wong, and P. Khatri, "A comprehensive time-course-based multicohort analysis of sepsis and sterile inflammation reveals a robust diagnostic gene set," *Science translational medicine*, vol. 7, no. 287, pp. 287ra71–287ra71, 2015.
- [6] T. E. Sweeney, H. R. Wong, and P. Khatri, "Robust classification of bacterial and viral infections via integrated host gene expression diagnostics," *Science translational medicine*, vol. 8, no. 346, pp. 346ra91–346ra91, 2016.
- [7] T. E. Sweeney, T. M. Perumal, R. Henao, M. Nichols, J. A. Howrylak, A. M. Choi, J. F. Bermejo-Martin, R. Almansa, E. Tamayo, E. E. Davenport *et al.*, "A community approach to mortality prediction in sepsis via gene expression analysis," *Nature communications*, vol. 9, no. 1, p. 694, 2018.
- [8] T. Barrett, S. E. Wilhite, P. Ledoux, C. Evangelista, I. F. Kim, M. Tomashevsky, K. A. Marshall, K. H. Phillippy, P. M. Sherman, M. Holko, A. Yefanov, H. Lee, N. Zhang, C. L. Robertson, N. Serova, S. Davis, and A. Soboleva, "NCBI GEO: archive for functional genomics data sets—update," *Nucleic Acids Research*, vol. 41, no. D1, pp. D991–D995, 11 2012. [Online]. Available: https://doi.org/10.1093/nar/gks1193
- [9] A. Brazma, H. Parkinson, U. Sarkans, M. Shojatalab, J. Vilo, N. Abeygunawardena, E. Holloway, M. Kapushesky, P. Kemmeren, G. G. Lara *et al.*, "ArrayExpress—a public repository for microarray gene expression data at the EBI," *Nucleic acids research*, vol. 31, no. 1, pp. 68–71, 2003.
- [10] M. B. Mayhew, L. Buturovic, R. Luethy, U. Midic, A. R. Moore, J. A. Roque, B. D. Shaller, T. Asuni, D. Rawling, M. Remmel *et al.*, "A generalizable 29-mrna neural-network classifier for acute bacterial and viral infections," *Nature communications*, vol. 11, no. 1, p. 1177, 2020.

- [11] scikit learn, "sklearn.model_selection.groupkfold," https://scikit-learn.org/stable/modules/ generated/sklearn.model_selection.GroupKFold.html, 2024, accessed: 2024-04-04.
- [12] D. J. Hand and R. J. Till, "A simple generalisation of the area under the roc curve for multiple class classification problems," *Machine learning*, vol. 45, pp. 171–186, 2001.
- [13] B. J. Parker, S. Günter, and J. Bedo, "Stratification bias in low signal microarray studies," BMC bioinformatics, vol. 8, pp. 1–16, 2007.
- [14] A. Airola, T. Pahikkala, W. Waegeman, B. De Baets, and T. Salakoski, "A comparison of auc estimators in small-sample studies," in *Machine learning in systems biology*. PMLR, 2009, pp. 3–13.
- [15] T. Akiba, S. Sano, T. Yanase, T. Ohta, and M. Koyama, "Optuna: A next-generation hyperparameter optimization framework," in *Proceedings of the 25th ACM SIGKDD international conference on knowledge discovery & data mining*, 2019, pp. 2623–2631.
- [16] M. B. Mayhew, U. Midic, K. Choi, P. Khatri, L. Buturovic, and T. E. Sweeney, "Towards equitable patient subgroup performance by gene-expression-based diagnostic classifiers of acute infection," *medRxiv*, pp. 2022–04, 2022.
- [17] United States Food and Drug Administration, "Cfr code of federal regulations title 21," https: //www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm?CFRPart=820, 2023, accessed: 2024-05-15.
- [18] S. R. Hayden and M. D. Brown, "Likelihood ratio: a powerful tool for incorporating the results of a diagnostic test into clinical decision making," *Annals of emergency medicine*, vol. 33, no. 5, pp. 575–580, 1999.
- [19] United States Food and Drug Administration, "SeptiCyte RAPID 510(k) Decision Summary," https://www.accessdata.fda.gov/cdrh_docs/reviews/K203748.pdf, 2021, accessed: 2024-05-12.
- [20] R. Luethy and L. Buturovic, "Decision Threshold Optimization for Diagnostic Tests using a Genetic Algorithm," https://inflammatix.com/decision-threshold-optimization-blog, 2023, accessed: 2024-05-01.
- [21] F.-A. Fortin, F.-M. De Rainville, M.-A. G. Gardner, M. Parizeau, and C. Gagné, "Deap: Evolutionary algorithms made easy," *The Journal of Machine Learning Research*, vol. 13, no. 1, pp. 2171–2175, 2012.
- [22] J. Ducharme, W. H. Self, T. M. Osborn, N. A. Ledeboer, J. Romanowsky, T. E. Sweeney, O. Liesenfeld, and R. E. Rothman, "A multi-mrna host-response molecular blood test for the diagnosis and prognosis of acute infections and sepsis: proceedings from a clinical advisory panel," *Journal of personalized medicine*, vol. 10, no. 4, p. 266, 2020.
- [23] Y. D. He, E. M. Wohlford, F. Uhle, L. Buturovic, O. Liesenfeld, and T. E. Sweeney, "The optimization and biological significance of a 29-host-immune-mrna panel for the diagnosis of acute infections and sepsis," *Journal of Personalized Medicine*, vol. 11, no. 8, p. 735, 2021.

8 Technical Appendices

8.1 Software architecture

The Myrna[™] device runs Windows OS. The software architecture, including classifiers, is shown in Fig. 4.

8.2 Datasets for training and validation

Development data are summarized in Tables 3 and 4. Healthy controls are not listed in the tables because the detailed demographic information was not available.



Figure 4: **The system software overview.** BVN: Bacterial/Viral/Noninfected; SEV: Severity; SP: Signal Processing; LAMP: Loop-Mediated Isothermal Amplification.

Table 3: Development datasets for diagnostic classifier (bacterial, viral, noninfected). ICU:
Intensive Care Unit. IQR: Interquartile Range. Unk.: Unknown. COPD: Chronic Obstructive Pul-
monary Disease. CAP: Community-Acquired Pneumonia. SIRS: Systemic Inflammatory Response
Syndrome.

Juniy	Autioi	Description	18	Age	winc	Tiationin	Country	Ducteriui	• • • •	Ron-Infected
		Patients hospitalized with COPD								10.00
EMEXP3589	Almansa	exacerbation	27	Unk.	16 (59)	Agilent	Spain	4 (15)	5 (19)	18(67)
		Surgical patients with sepsis								
EMTAB1548	Almansa	(EXPRESS)	155	72.0 (IQR 61.0-78.0)	95 (61)	Agilent	Spain	82 (53)	0(0)	73 (47)
EMTAB3162	van den Ham	Patients with dengue	21	20.0 (IQR 17.0-28.0)	11 (52)	Affymetrix	Indonesia	0 (0)	21 (100)	0 (0)
		Sepsis due to faecal								
EMTAB5273	Burnham	peritonitis or pneumonia	129	66.0 (IQR 51.0-76.0)	59 (46)	Illumina	United Kingdom	119 (92)	0(0)	10 (8)
		Sepsis due to faecal								
EMTAB5274	Burnham	peritonitis or pneumonia	108	71.0 (IQR 62.75-77.0)	69 (64)	Illumina	United Kingdom	108 (100)	0 (0)	0 (0)
EMTAB5638	Almansa	Ventilator-associated pneumonia in ICU	17	Unk.	0 (0)	Agilent	Spain	0 (0)	0 (0)	17 (100)
00540004		Severe H1N1/09 influenza or	02	1.25 (100) 0.20 (1.0)	22.00		N 5 1 K 1	10 (20)	11.(15)	22.00
GSE42026	Herberg	bacterial infection	92	1.25 (IQR 0.38-4.0)	33 (36)	Illumina	United Kingdom	18 (20)	41 (45)	33 (56)
GSE42834	Bloom	Bacterial pneumonia or sarcoidosis	200	Unk.	32 (16)	Illumina	United Kingdom	14 (7)	0(0)	186 (93)
GSE47655	Stone	Severe anaphylaxis in ED	12	Unk.	0(0)	Affymetrix	Australia	0(0)	0(0)	12 (100)
GSE51808	Kwissa	Patients with dengue	37	Unk.	0(0)	Affymetrix	Thailand	0(0)	28 (76)	9 (24)
GSE57065	Cazalis	Septic shock	53	62.0 (IQR 54.25-76.0)	19 (36)	Affymetrix	France	28 (53)	0(0)	25 (47)
GSE57183	Shenoi	Febrile children with SJIA	14	3.6 (IQR 3.3-7.3)	5 (36)	Illumina	United States of America	0(0)	0(0)	14 (100)
GSE25504	Smith	Neonatal sepsis	39	0.0 (IQR 0.0-0.0)	19 (49)	Affymetrix	United Kingdom	29 (74)	4 (10)	6 (15)
GSE60244	Suarez	Lower respiratory tract infections	133	63.0 (IQR 50.0-77.0)	37 (28)	Illumina	United States of America	22(17)	71 (53)	40 (50)
		innuenza patients or varying								
GSE61821	Hoang	seventy	48	40.0 (IQR 19.75-51.0)	24 (50)	Illumina	Vietnam	0(0)	48 (100)	0 (0)
		Pediatric patients with Kawasaki								
GSE63881	Hoang	disease	171	2.75 (IQR 1.42-4.25)	102 (60)	Illumina	United States of America	0(0)	0(0)	171 (100)
CEEGAASG	Mahaing	Febrile infants with bacterial	210	0.1 (COR 0.06 0.12)	106 (49)	111	United States of America	80 (41)	111 (51)	10 (0)
GSE04450	Manajan	or viral infection in ED Sucpacted but negative for	219	0.1 (IQR 0.06-0.13)	106 (48)	IIIumina	United States of America	89 (41)	111 (51)	19 (9)
CSE65692	Saichuna	CAP	75	50.0 (108.48.0.67.0)	22 (28)	Affamatria	Natharlands	0.00	0.00	75 (100)
G5E05082	Sciciuna	LAP Padiataia ICU nationta with	/5	.59.0 (IQK 48.0-07.0)	22 (28)	Allymetrix	inetheriands	0(0)	0(0)	/5 (100)
CSE66000	Swaanay	rediatric ICU patients with	107	2.45 (TOP 1.0.5.99)	04 (48)	Affamatria	United States of America	100 (55)	11.60	77 (30)
GSE00099	Sweeney	sepsis, SIKS	19/	2.45 (IQK 1.0-5.88)	52 (52)	Anymetrix	United States of America	109 (55)	11(0)	21 (39)
GSE67059	Heinonen	Children with rhinovirus	101	0.83 (IQR 0.3-1.29)	53 (52)	Illumina	United States of America	0(0)	80 (79)	21 (21)
L D C HORE	M 101 - 1	Trauma patients with or	272	22.0.000.21.0.11.0	00.000			47 (12)	0.00	225 (07)
glueBuffyHCSS	Multiple authors	without infection	372	32.0 (IQR 24.0-41.0)	98 (26)	Affymetrix	United States of America	47(13)	0(0)	325 (87)
		Outpatients with acute respiratory								
GSE68310	Zhai	viral infections	347	20.96 (IQR 20.09-22.76)	50 (14)	Illumina	United States of America	0(0)	104 (30)	243 (70)
		Sepsis, many cases from								
GSE69528	Khaenam	burkholderia	138	Unk.	0(0)	Illumina	Thailand	83 (60)	0(0)	55 (40)
		Pediatric patients with bacteriai			0.000					
GSE/2810	Herberg	or viral infection	15	1.83 (IQR 0.88-3.29)	8 (53)	Illumina	United Kingdom	5 (33)	10(67)	0 (0)
		Pediatric patients with bacterial								
GSE/3461	Wright	or viral infection or other inflammatory diseases	363	2.79 (IQR 0.92-8.81)	165 (45)	Illumina	United Kingdom	52 (14)	94 (26)	217 (60)
		Infants with respiratory synctial								
GSE77087	de Steenhuijsen Piters	virus	59	0.45 (IQR 0.14-0.69)	25 (42)	Illumina	United States of America	0 (0)	41 (69)	18 (31)
GSE77791	Plassais	Severe burn shock	30	48.0 (IQR 40.25-55.0)	21 (70)	Affymetrix	France	0 (0)	0 (0)	30 (100)
CEEROOSO	Terre	Moderate and severe influenza	20	64 5 (IOD 48 5 74 25)	12 (22)	A millions	6	0.00	24 (62)	15 (29)
GSE82050	Tang	infection	- 39	64.5 (IQR 48.5-74.25)	13 (33)	Agiient	Germany	0(0)	24 (02)	15 (38)
0000100000	B 11 B 1	infants with respiratory syncular	74	0.05 (100 0.15 0.10)	20 (52)		11 1 10 1 11 11	0.00	(2.00)	10.00
GSE103842	Rodriguez-Fernandez	Virus	74	0.25 (IQR 0.17-0.44)	39 (53)	Illumina	United States of America	0(0)	62 (84)	12(16)
GSEI11368	Dunning	Innuenza HINI and B	33	38.0 (IQR 29.0-49.0)	15 (45)	Illumina	United Kingdom	0(0)	33 (100)	0(0)
									0.00	
000000000000000000000000000000000000000	D 11	Sepsis, many cases from		54.0 (100) 40.0 (1.0)	26 (17)	Y11 -	ann 11 1	45 (02)		10 (10)
GSE13015gpl6102	Pankla	Sepsis, many cases from burkholderia	55	54.0 (IQR 48.0-61.0)	26 (47)	Illumina	Thailand	45 (82)	0(0)	10 (18)
GSE13015gpl6102	Pankla	Sepsis, many cases from burkholderia Sepsis, many cases from burkholderia	55	54.0 (IQR 48.0-61.0)	26 (47)	Illumina	Thailand	45 (82)	0 (0)	10 (18)
GSE13015gpl6102 GSE13015gpl6947	Pankla Pankla	Sepsis, many cases from burkholderia Sepsis, many cases from burkholderia	55 20	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5)	26 (47) 6 (30)	Illumina	Thailand	45 (82) 15 (75)	0(0)	10 (18) 5 (25)
GSE13015gpl6102 GSE13015gpl6947 GSE21802	Pankla Pankla Bermejo-Martin	Sepsis, many cases from burkholderia Sepsis, many cases from burkholderia Pandemic H1N1 in ICU	55 20 14	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk.	26 (47) 6 (30) 0 (0)	Illumina Illumina Illumina	Thailand Thailand Unk.	45 (82) 15 (75) 0 (0)	0 (0) 10 (71)	10 (18) 5 (25) 4 (28)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22008	Pankla Pankla Bermejo-Martin	Sepsis, many cases from burkholderia Sepsis, many cases from burkholderia Pandemic HINI in ICU Bacterial infection or other informatics conditione	55 20 14	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk.	26 (47) 6 (30) 0 (0) 59 (22)	Illumina Illumina Illumina	Thailand Thailand Unk.	45 (82) 15 (75) 0 (0) 52 (10)	0 (0)	10 (18) 5 (25) 4 (28) 222 (81)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE227131	Pankla Pankla Bermejo-Martin Berry Bardal	Sepas, mary cases from burkholderia Sepais, many cases from burkholderia Pandemic HINI in ICU Bacterial infection or other inflammatory conditions	55 20 14 274	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75.44.75)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71)	Illumina Illumina Illumina Illumina	Thailand Thailand Unk. Unk.	45 (82) 15 (75) 0 (0) 52 (19) 0 (0)	0 (0) 10 (71) 0 (0) 7 (50)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE22730	Pankla Pankla Bermejo-Martin Bertal Sutharland	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Pandemic HINI in ICU Bacterial Infection or other inflammatory condition or other Severe HINI Severe HINI	55 20 14 274 14	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk.	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0)	Illumina Illumina Illumina Affymetrix	Thailand Thailand Unk. Unk. Norway Australia	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (49)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE27131 GSE28750 GSE28750	Pankla Pankla Bermejo-Martin Berry Berdal Sutherland	Sepas, mary cases from burkholderia nonnari Separation and the second second Pandemic HINI in GCU Bacterial Infection or other inflammatory conditions Secret HINI Sepais or post-angleal SIRS	55 20 14 274 14 21	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. Unk.	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0)	Illumina Illumina Illumina Affymetrix Affymetrix Illuming	Thailand Thailand Unk. Unk. Norway Australia Unk	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE27131 GSE28750 GSE28750 GSE28991	Pankla Pankla Bermejo-Martin Berdal Berdal Sutherland Naim	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Pandemic HINI in ICU Bacterial infection or other inflammatory conditions Severe HINI Sepais or poss-surgical SIRS Pattents with saspected dengue	55 20 14 274 14 21 11	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. Unk.	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina	Thailand Thailand Unk. Unk. Norway Australia Unk.	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE2708 GSE27131 GSE28750 GSE287991 GSE28991	Pankla Pankla Bermejo-Martin Berry Berdal Sutherland Naim	Sepas, mary cases from burkholder, from Sepais, mary cases from Pendemic HINI in ICU Bacierial infection or other inflammatory conditions Severe HINI Sepais or post-surgical SIRS Patients with suspected dengue Febrie patients with viral	55 20 14 274 14 21 11 80	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. 25.0 (IQR 22.0 40.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina	Thailand Thailand Unk. Unk. Norway Australia Unk.	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 0 (0)
GSE13015gpl6402 GSE13015gpl6497 GSE21802 GSE22098 GSE27131 GSE28750 GSE28991 GSE29385	Pankla Pankla Bermejo-Martin Berry Berdal Sutherland Naim Naim	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Pandemic HINI in ICU Bacterial infection or other inflammatory conditions Severe HINI Sepais or poss-surgical SIRS Pattens with suspected deigue Febrile patients with viral Bacteria, marky influenza	55 20 14 274 14 21 11 80	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. 0.00	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina	Thailand Thailand Unk. Unk. Norway Australia Unk. Unk.	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 0 (0)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE27131 GSE28750 GSE28991 GSE28991 GSE29385	Pankla Pankla Bermejo-Martin Berry Berdal Sutherland Naim Naim	Sepas, mary cases from burkholderia Sepais, mary cases from Pendemic HINS in a CU Basieral Infection or other infammatory conditions Severe HINI Sepais or post-surgical SIRS Patients with surgected dengue Febrile patients with viral infection, mostly influenza Pediatric patients with viral	55 20 14 274 14 21 11 80 81	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. Unk. 25.0 (IQR 22.0-40.0) 55.0 (IQR 12.0 U.C)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 24 (42)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina	Thailand Thailand Unk. Unk. Norway Australia Unk. Unk.	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0) 50 (72)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 0 (0) 22 (23)
GSE13015gpl6102 GSE13015gpl6047 GSE21802 GSE21802 GSE21802 GSE27131 GSE28750 GSE28750 GSE29385 GSE30119	Pankla Pankla Bernejo-Martin Berty Berdal Sutherland Naim Naim Banchereau	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Pandemic HINI in ICU Bacterial infection or other inflammatory conditions Severe HINI Sepais or poss-surgical SIRS Patients with suspected deingue Febrile patients with viral Infection, mostly influenza Pediatric patients with community-sequired bacterial infection marking informa-	55 20 14 274 14 21 11 80 81	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. 25.0 (IQR 22.0-40.0) 6.5 (IQR 1.92-11.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 34 (42)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina	Thailand Thailand Unk. Unk. Norway Australia Unk. Unk. Unk.	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0) 59 (73)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 0 (0) 22 (27)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE2750 GSE28750 GSE28750 GSE29385 GSE30119 CSE23757	Pankla Pankla Bernejo-Martin Berrdal Sutherland Naim Naim Banchereau	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Bacterial network of the second Bacterial network of the second Severe HINI Sepais or post-surgical SIRS Patients with surgeted dengue Febrile patients with viral infection, mostly influenza Pedatire patients with community-acquired bacterial infection Patients undergoing mechanical ventilation	55 20 14 274 14 21 11 80 81 44	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.754.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 6.5 (0QR 1.6 6.5.2)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 12 (20)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina	Thailand Thailand Unk. Norway Australia Unk. United States of America	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0) 59 (73) 0 (0)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 0 (0) 22 (27) 44 (100)
GSE13015gpl6102 GSE13015gpl6047 GSE21802 GSE21802 GSE22008 GSE22703 GSE28750 GSE2891 GSE29385 GSE30119 GSE32707	Pankla Pankla Bernejo-Martin Berdy Berdal Sutherland Naim Naim Banchereau Dolinay	Sepas, mary cases from burkholderia or the statistical of the statistical of the statistical particular statistical of the statistical particular statistical of the statistical statistical statistical statistical statistical statistical statistical statistical statistical particular statistical statistical statistical particular statistical	55 20 14 274 14 21 11 80 81 44	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9.2-11.0) 56.0 (1QR 45.0-59.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30)	Illumina Illumina Illumina Aflymetrix Aflymetrix Illumina Illumina Illumina	Thailand Thailand Unk. Unk Norway Australia Unk. Unk. United States of America United States of America	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0) 59 (73) 0 (0)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE22098 GSE28750 GSE28750 GSE29991 GSE29991 GSE30119 GSE30119	Pankla Pankla Bermejo-Martin Berty Berdal Sutherland Naim Naim Banchereau Dolinay Parnell	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Pandemie IINI in CI Pandemie IINI in CI Pandemie IINI in CI Sepais or post-surgical SIRS Patients with surged SIRS Patients with surged SIRS Patients with original refection, mostly influenza Pediatric patients with community-acquired bacterial infection Patients undergoing mechanical ventilation with sepais or SIRS or neither Bacterial or milleraz A	55 20 14 274 14 21 11 80 81 44 54	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9-21.10.0) 56.0 (1QR 45.0-59.0) 50.0 (1QR 45.0-59.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina	Thailand Thailand Unk. Unk. Unk. Unk. Unk. United States of America United States of America Australia	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 8 (15)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100) 30 (56)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22103 GSE22098 GSE27131 GSE2891 GSE2891 GSE29385 GSE30119 GSE32707 GSE40012 GSE40015	Pankla Pankla Bernejo-Martin Berty Bertdal Sutherland Naim Banchereau Dolinay Panell Nawen	Sepas, mary cases from burkhoders, many cases Sepais, mary cases Pandemie HIN in fCU Bacieral infection or other inflammatory conditions Severe HINI Septier and Severe HINI Severe HINI S	55 20 14 274 14 21 11 80 81 44 54 123	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-57.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina	Thailand Thailand Unk. Unk. Unk. Unk. Unk. Unk. United States of America Australia Viennom	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 8 (15) 123 (100)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100) 30 (56) 0 (0)
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22088 GSE22078 GSE22078 GSE28750 GSE28750 GSE28750 GSE2985 GSE32707 GSE32707 GSE32707 GSE40012 GSE40015	Pankla Pankla Bermejo-Martin Berdy Berdal Sutherland Naim Naim Banchereau Dolinay Panell Nguyen	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Pendemic HINI in ICU Pendemic HINI in ICU Sepais or post-surgical SIRS Patients with surgected dengue Febrie patients with viral infection, mostly influenza Pedatire jutants with community-acquired bacterial infection Patients undergoing mechanical ventilation with explois or SIRS or neither Bacterial or SIRS or neither Bacterial ventilemaza A penemonia or SIRS penemonia or SIRS penemonia or SIRS	55 20 14 274 14 21 11 80 81 44 54 123	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 45.6-70.0) 12.0 (1QR 10.0-14.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina	Thailand Thailand Unk. Unk. Norway Australia Unk. United States of America United States of America United States of America	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 8 (15) 123 (100)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100) 30 (56) 0 (0)
GSE13015gpl6102 GSE13015gpl6047 GSE21802 GSE21802 GSE22098 GSE22098 GSE2750 GSE28931 GSE29385 GSE30119 GSE32707 GSE40102 GSE40105 GSE40165	Pankla Pankla Bernyio-Martin Berdal Sutherland Naim Banchereau Dolinay Parnell Nguyen Hu	Sepas, mary cases from burkholder, Sepais, mary cases from Pendemie HIN in CU Bacieral infection or other inflammatory conditions Severe HINI Sepsis or post-surgical SIRS Patients with suspected dengue Febrie patients with viral Pediatic patients with viral Pediatic patients with community-sequence bacterial infection Patients undergoing mechanical ventilation where the sepsis or SIRS or neither Bacterial or influenza A pneumonia or SIRS Febrile children with dengua Febrile Children with dengua	55 20 14 274 14 21 11 80 81 44 54 123 52	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 15.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 55.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-57.0) 12.0 (1QR 11.0-114.0) 0.92 (1QR 0.33.1 6)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina	Thailand Thailand Unk. Unk. Unk. Unk. Unk. Unk. United States of America Australia Vietnam United States of America	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0) 8 (15)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 8 (15) 123 (100) 22 (42)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100) 30 (56) 0 (0) 22 (42)
CSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22008 GSE27131 GSE22988 GSE20119 GSE32077 GSE32077 GSE32077 GSE40012 GSE40015 GSE40165	Pankla Pankla Bernejo-Martin Berdal Subterland Naim Naim Danchereau Dolinay Pamell Nguyen Hu Lill	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Beitreila Infection future Beitreila Infection official Severe HINI Sepais or post-surgical SIRS Patients with surgeted dengue Febrie patients with viral infection, mostly influenza Pedatric patients with viral infection mostly influenza Pedatric patients with community-acquired bacterial infection Patients undergoing mechanical ventilation with sepsis or SIRS or neither Bacterial or influenza A pneumonia or SIRS Pedric enlidere with dengue reference information and the dengue or viral infection	55 20 14 274 14 21 11 80 81 44 54 123 52 21	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-57.0) 59.0 (1QR 45.0-57.0) 59.0 (1QR 45.0-70.0) 57.0 (1QR 33.0-70.5) 57.0 (1QR 33.0-70.5)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Affymetrix	Thailand Thailand Unk. Norway Australia Unk. United States of America Australia Vietnam United States of America Linited States of America	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0) 8 (15) 21 (100)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100) 30 (56) 0 (0) 22 (42) 0 (0)
CSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE21802 GSE22908 GSE29985 GSE29991 GSE29191 GSE29191 GSE29191 GSE29191 GSE29191 GSE29102 GSE30119 GSE30119 GSE40012 GSE40012 GSE40356	Pankla Bernejo-Martin Berry Berdal Sutherland Naim Banchereau Dolinay Pannell Nguyen Hu Lall	Sepas, mary cases from burkholderia Sepais, mary cases from Pendemie HIN in a CU Bacieral infaction or other infammatory conditions Severe HIN Sepaie and Severe HIN Severe HIN	55 20 14 274 14 21 11 80 81 44 54 123 52 21	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 56.0 (1QR 45.0-70.0) 12.0 (1QR 10.0-14.0) 9.0 (1QR 45.5-67.0) 12.0 (1QR 13.16) 57.0 (1QR 33.0-70.5)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina	Thailand Thailand Unk. Unk. Unk. United States of America Vietnam Vietnam United States of America Estonia	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0) 8 (15) 21 (100)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 0 (0) 0 (0) 20 (0) 22 (42) 0 (0)	$\begin{array}{c} 10 \ (18) \\ \hline \\ 5 \ (25) \\ 4 \ (28) \\ \hline \\ 222 \ (81) \\ 7 \ (50) \\ 11 \ (52) \\ 0 \ (0) \\ \hline \\ 22 \ (27) \\ \hline \\ 44 \ (100) \\ \hline \\ 30 \ (56) \\ 0 \ (0) \\ \hline \\ 22 \ (42) \\ 0 \ (0) \\ \hline \end{array}$
GSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE21802 GSE22008 GSE22008 GSE27530 GSE29385 GSE2071 GSE2077 GSE2077 GSE2077 GSE30119 GSE30105 GSE4012 GSE40165 GSE40165 GSE40586 INF-02	Pankla Pankla Bermejo-Martin Berdal SutberInd Naim Naim Banchereau Dolinay Parnell Nguyen Hu Läll Liesenfeld	Sepas. mary cases from burkholderia Punderia (Sepandi Sepandi Sepandi Sepandi Pundemie HINI in ICU Bacterial Infection or other inflammatory conditions Server HINI Sepandi Sepandi Sepandi Sepandi Sepandi Sepandi Sepandi Sepandi Sepandi Infection, mosty influenza Pechatric patients with viral infection, mosty influenza Pechatric patients with community-sequired bacterial infection Patients undergoing mechanical ventilation with sepais or SIRS or neither peneumonia or SIRS Pebrile children with bacterial or viral infection Community-sequired bacterial meringitis Pelorine children with bacterial or viral infection Patients undergoing in the sepandi Sepandi	55 20 14 274 41 21 11 80 81 44 54 54 52 21 59	54.0 (IQR 48.0-61.0) 49.0 (IQR 43.5-59.5) Unk. 16.0 (IQR 11.0-26.0) 39.0 (IQR 35.75-44.75) Unk. Unk. 25.0 (IQR 22.0-04.0) 6.5 (IQR 1.92-11.0) 56.0 (IQR 45.0-59.0) 59.0 (IQR 45.0-59.0) 0.92 (IQR 0.33-1.6) 57.0 (IQR 30.3-70.5) 9.0 (IQR 20.5-50.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Mumina NanoStrine®	Thailand Thailand Unk. Nerway Australia Unk. Unk. United States of America Australia Vietnam United States of America Lunited States of America Estonia United States of America Estonia	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0) 16 (30) 0 (0) 21 (100) 24 (41)	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 7 (50) 11 (100) 80 (100) 0 (0) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0) 7 (12)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 22 (27) 44 (100) 30 (56) 0 (0) 22 (42) 0 (0) 28 (47)
CSE13015gpl6102 GSE13015gpl6047 GSE21802 GSE21802 GSE22098 GSE22098 GSE2985 GSE2985 GSE23707 GSE23707 GSE23707 GSE23707 GSE30119 GSE3017 GSE40012 GSE40016 GSE40396 INF-02	Pankla Bernejo-Martin Berry Berdal Sutherland Naim Banchereau Dolinay Pamell Nguyen Hu Lili Liesenfeld	Sepas, may cases from burkholderia Sepais, may cases from Pendemic HINE in a CU Bacieral Infection or other infammatory conditions Severe HINI Sepais or post-surgical SIRS Patients with surgected dengue Febrile patients with oral infection, mostly influenza Pediatric patients with oral patients with community-acquired bacterial infection material infection patients with community-acquired bacterial or influenza A perturbation or SIRS Febrile children with dengue Febrile children with a the febrile KCU patients with at	55 20 14 274 14 21 11 80 81 44 54 123 52 21 59	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9-211.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-70.0) 12.0 (1QR 10.0-14.0) 0.92 (1QR 0.33.1-6) 57.0 (1QR 35.0-70.5) 39.0 (1QR 29.0-56.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina Illumina NanoString [®]	Thailand Thailand Unk. Unk. Unk. United States of America United States of America Australia Vietnam United States of America Australia United States of America Estonia	$\begin{array}{c} 45 \ (82) \\ 15 \ (75) \\ 0 \ (0) \\ 52 \ (19) \\ 0 \ (0) \\ 52 \ (19) \\ 0 \ (0) \\ 0 \ (0) \\ 59 \ (73) \\ 0 \ (0) \\ 59 \ (73) \\ 0 \ (0) \\ 16 \ (30) \\ 0 \ (0) \\ 8 \ (15) \\ 21 \ (100) \\ 24 \ (41) \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 123 (100) 22 (42) 0 (0) 7 (12)	$\begin{array}{c} 10 \ (18) \\ \hline 10 \ (18) \\ \hline 5 \ (25) \\ \hline 4 \ (28) \\ \hline 7 \ (50) \\ \hline 7 \ (50) \\ \hline 0 \ (0) \\ \hline 0 \ (0) \\ \hline 22 \ (27) \\ \hline 44 \ (100) \\ \hline 30 \ (56) \\ \hline 0 \ (0) \\ \hline 22 \ (42) \\ \hline 0 \ (0) \\ \hline 28 \ (47) \\ \hline \end{array}$
GSE13015gpl6102 GSE13015gpl6027 GSE21802 GSE21802 GSE22098 GSE22098 GSE2985 GSE2985 GSE30119 GSE30119 GSE30119 GSE30119 GSE4015 GSE40165 GSE40165 GSE40186 GSE40185 GSE40185 GSE40185 GSE40185 GSE40185 GSE40185 GSE40186 GSE40185	Pankla Pankla Bernejo-Martin Berry Berdal Sutherland Naim Naim Banchereau Dolinay Panell Nguyen Hu Lill Liesenfeld Rogers	Sepas, may cases from burkhoderia Brand and the sense of the sense Pandemic HIN in RCU Bacterial Infection or other inflammatory conditions Severe HINI Septement Sense and Sense Patterns with suspected dengue infection, mostly inflamena bacterial infection Patients undergoing mechanical ventilation with sepsisor SIRS or neither Bacterial or influenza Bacterial or influenza Bacterial or influenza Febric children with bacterial or viral Infection Community-sequired bacterial meningitis Patients with subjected inflamenta Community-sequired bacterial meningitis Patients with subjected inflamenta Community-sequired bacterial meningitis Patients with subjected inflamentary of the CULP patients with a fease Batteria for far actar ensitient or distance sendomic Information of subjects and the sendomic of the second Batteria for far actar ensitient or distance sendomic Batteria for far ensitient ensitient or distance sendomic Batteria for far actar ensitient or distance sendomic Batteria for action far ensitient or distance sendomic Batteria for far ensit ensitient or	55 20 14 274 14 21 11 80 81 44 54 123 52 21 59 13	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 12.0-11.0) 50.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-79.0) 59.0 (1QR 0.33-1.6) 57.0 (1QR 24.0-77.0) 57.0 (1QR 24.0-77.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77)	Illumina Illumina Illumina Alfymetrix Alfymetrix Illumina Illumina Illumina Illumina Illumina NamoString® NanoString®	Thailand Thailand Unk. Norway Australia Unk. Unk. United States of America Australia Vietnam United States of America Estonia United States of America Estonia	$\begin{array}{c} 45 (82) \\ 15 (75) \\ 0 (0) \\ 52 (19) \\ 0 (0) \\ 52 (19) \\ 0 (0) \\ 0 (0) \\ 0 (0) \\ 0 (0) \\ 59 (73) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 21 (100) \\ 21 (100) \\ 24 (41) \\ 6 (46) \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 20 (0) 2123 (100) 22 (42) 0 (0) 7 (12) 1 (8)	10(18) 5(25) 4(28) 7(50) 11(52) 0(0) 0(0) 22(27) 44(100) 30(56) 0(0) 22(42) 0(0) 22(42) 0(0) 24(47) 6(46)
CSE13015gpl6102 GSE13015gpl6047 GSE21802 GSE21802 GSE22098 GSE29385 GSE29385 GSE20707 GSE20707 GSE20707 GSE20707 GSE20707 GSE40012 GSE40016 GSE40056 INF-02 INF-01	Pankla Penkla Berrncjo-Martin Berry Berdal Sutherland Naim Naim Banchereau Dolinay Parmell Nguyen Hu Luiseenfeld Rogers	Sepas, may cases from burkholderia Sepais, may cases from burkholderia Bacietai Interior or other infammatory conditions Severe HINI Sepais or post-augical SIRS Patients with suspected dengue Febrile patients with viral infection, mostly influenza Pedatre patients with community-acquired bacterial infection with sepais or SIRS or nother Patients using or SIRS Pethele children with dengue Febrile children dengue Febrile	55 20 14 274 14 21 11 80 81 44 54 52 21 59 13	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9-21.1.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 46.5-67.0) 12.0 (1QR 10.0-14.0) 0.92 (1QR 0.33-1.6) 57.0 (1QR 33.0-70.5) 39.0 (1QR 29.0-56.0) 57.0 (1QR 23.0-72.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77)	Illumina Illumina Illumina Illumina Alfymetrix Alfymetrix Illumina Illumina Illumina Illumina Illumina Alfymetrix NanoString®	Thailand Thailand Unk Norway Australia Unk Unk United States of America Australia Vietnam United States of America Estonia USAGreece USA	$\begin{array}{c} 45 \ (82) \\ 15 \ (75) \\ 0 \ (0) \\ 52 \ (19) \\ 0 \ (0) \\ 52 \ (19) \\ 0 \ (0) \\ 0 \ (0) \\ 59 \ (73) \\ 0 \ (0) \\ 59 \ (73) \\ 0 \ (0) \\ 16 \ (30) \\ 0 \ (0) \\ 8 \ (15) \\ 21 \ (100) \\ 24 \ (41) \\ 6 \ (46) \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 80 (100) 0 (0) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0) 7 (12) 1 (8)	$\begin{array}{c} 10(18)\\ 5(25)\\ 4(28)\\ 222(81)\\ 7(50)\\ 11(52)\\ 0(0)\\ 22(27)\\ 44(100)\\ 30(56)\\ 0(0)\\ 22(42)\\ 0(0)\\ 22(47)\\ 6(46)\\ \end{array}$
GSE13015gp16102 GSE13015gp16947 GSE21802 GSE22998 GSE27131 GSE2985 GSE30119 GSE32907 GSE30119 GSE3207 GSE40012 GSE40012 GSE40012 GSE40956 GSE40957	Pankla Pankla Bernyi Berdal Sutherland Naim Bancherau Dolinay Parnell Nguyen Hu Lill Liesenfeld I isenfeld	Sepas, mary cases from burkholder, Pendemie HINI in CU Pandemie HINI in CU Bacieral infection or other inflammatory conditions Severe HINI Seysis or post-surgical SIRS Patients with suspected dengue Febrile patients with train Petiatric patients with train Petiatric patients with community-sequence bacterial infection Patients undergoing mechanical ventilation with sepsis or SIRS or neither Bacterial or influenza A patients undergoing mechanical ventilation with sepsis or SIRS or neither Bacterial or influenza A patients undergoing mechanical ventilation vital infection Community-sequence bacterial meningitis Patients with suspected inflections/ Patients with a flections/ Patients with a flection flections/ Patients with a flection flections/ Patients with a flection flections/ Patients with a flection flections/ Patients with flection flections/ Patients with flections/ Patients with flections/ Patients with flections/ Patients/	55 20 14 274 14 21 11 80 81 44 54 123 52 21 59 13	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (UQR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 50.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 57.0 (1QR 53.0-70.5) 39.0 (1QR 29.0-56.0) 57.0 (1QR 42.0-72.0) 56.0 (1QR 27.0-85.0) 57.0 (1QR 42.0-72.0) 56.0 (1QR 27.0-85.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62)	Illumina Illumina Illumina Aflymetrix Aflymetrix Illumina Illumina Illumina Illumina Illumina Illumina Illumina Marching [®] NanoString [®]	Thailand Thailand Unk. Unk. Unk. Unk. Unk. Unk. Unk. Unk.	$\begin{array}{c} 45 (82) \\ 15 (75) \\ 0 (0) \\ 52 (19) \\ 0 (0) \\ 52 (19) \\ 0 (0) \\ 0 (0) \\ 0 (0) \\ 0 (0) \\ 59 (73) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 21 (100) \\ 21 (100) \\ 24 (41) \\ 6 (46) \\ 5 (38) \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0) 7 (12) 1 (8) 8 (62)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 0 (0) 22 (27) 44 (100) 30 (56) 0 (0) 22 (42) 0 (0) 22 (42) 0 (0) 28 (47) 6 (46) 6 (46) 0 (0)
CSE13015gp16102 GSE13015gp16947 GSE21802 GSE22998 GSE29385 GSE29385 GSE29385 GSE29397 GSE29991 GSE29991 GSE29991 GSE30199 GSE30199 GSE40012 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE4	Pankla Bernejo-Martin Berry Berdal Sutherland Naim Banchereau Dolinay Panell Nguyen Hu Lill Liesenfeld Rogers Liesenfeld	Sepas, mary cases from burkholderia Sepais, mary cases from burkholderia Bactenine (INN) in COL Protection (INN) in COL Bactenine (INN) in COL Sepain (INN) in Col Sep	55 20 14 274 14 21 11 80 81 44 52 21 59 13	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9-211.0) 56.0 (1QR 45.0-59.0) 75.0 (1QR 45.67.0) 12.0 (1QR 10.0-14.0) 0.92 (1QR 0.33-1.6) 57.0 (1QR 35.07.05) 39.0 (1QR 42.0-72.0) 65.0 (1QR 42.0-72.0) 65.0 (1QR 37.0-85.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62)	Illumina Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Affymetrix NanoString® NanoString	Thailand Thailand Unk. Norway Australia Unk. United States of America United States of America Australia Vietnam United States of America Estonia USA/Greece USA	$\begin{array}{c} 45 \ (82) \\ 15 \ (75) \\ 0 \ (0) \\ 52 \ (19) \\ 0 \ (0) \\ 10 \ (48) \\ 0 \ (0) \\ 0 \ (0) \\ 59 \ (73) \\ 0 \ (0) \\ 16 \ (30) \\ 0 \ (0) \\ 8 \ (15) \\ 21 \ (100) \\ 24 \ (41) \\ 6 \ (46) \\ 5 \ (38) \end{array}$	3 (6) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0) 7 (12) 1 (8) 8 (62)	$\begin{array}{c} 10 \ (18) \\ \hline 10 \ (18) \\ \hline 5 \ (25) \\ \hline 4 \ (28) \\ \hline 222 \ (81) \\ \hline 7 \ (50) \\ \hline 11 \ (52) \\ \hline 0 \ (0) \\ \hline 0 \ (0) \\ \hline 22 \ (27) \\ \hline 44 \ (100) \\ \hline 30 \ (56) \\ \hline 0 \ (0) \\ \hline 22 \ (42) \\ \hline 0 \ (0) \\ \hline 28 \ (47) \\ \hline 6 \ (46) \\ \hline 0 \ (0) \\ \hline \end{array}$
CSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22898 GSE27181 GSE29895 GSE2985 GSE3019 GSE3077 GSE40012 GSE40012 GSE40012 GSE40056 INF-02 INF-1IS-01 INF-IIS-04	Pankla Pankla Bernejo-Martin Berry Berdal Sutherland Naim Naim Banchereau Dolinay Parnell Nguyen Hu Lift Liesenfeld Rogers Liesenfeld Giumpedlus-Boarbeutle	Sepas, mary cases from burkholderia Sepais, mary cases from Pendemie HINI in ICU Bacieral infection or other inflammatory conditions Severe HINI Sepais or post-surgical SIRS Patients with suspected dengue Febrile patients with viral bacterial inflamma pendemic severe HINI bacterial inflamma pendemic severe HINI bacterial inflamma pendemic severe severe bacterial inflamma pendemic severe severe bacterial inflamma pendemic severe bacterial inflamma pendemic severe bacterial inflamma pendemic severe bacterial inflamma pendemic severe bacterial inflamma pendemic severe bacterial inflamma pendemic severe pendemic severe	55 20 14 274 14 21 11 80 81 44 54 123 52 21 59 13 3	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (UQR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 55.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 57.0 (1QR 35.0-70.5) 39.0 (1QR 29.0-56.0) 57.0 (1QR 42.0-72.0) 65.0 (1QR 37.0-85.0) 72.0 (1QR 37.0-85.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62) 2 (67)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina Illumina Masostring® NanoString®	Thailand Thailand Unk. Unk. Unk. Unk. Unk. Unk. Unkets of America Australia Vietnam Vietnam United States of America Australia Vietnam United States of America Estonia USA/Greece USA	$\begin{array}{c} 45 (82) \\ 15 (75) \\ 0 (0) \\ 52 (19) \\ 0 (0) \\ 10 (48) \\ 0 (0) \\ 0 (0) \\ 59 (73) \\ 0 (0) \\ 59 (73) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 8 (15) \\ 21 (100) \\ 24 (41) \\ 6 (46) \\ 5 (38) \\ 3 (100) \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0) 7 (12) 1 (8) 8 (62) 0 (0)	10 (18) 5 (25) 4 (28) 222 (81) 7 (50) 11 (52) 0 (0) 2 (27) 44 (100) 30 (56) 0 (0) 22 (42) 0 (0) 22 (42) 0 (0) 28 (47) 6 (46) 0 (0) 0 (0)
GSE13015gpl6102 GSE10015gpl6947 GSE21802 GSE21802 GSE22908 GSE22908 GSE2985 GSE2985 GSE29810 GSE298102 GSE2985 GSE4012 GSE4015 GSE40165 GSE40165 INF-02 INF-18-03 INF-18-04	Pankla Pankla Bernejo-Martin Berry Berdal Sutherland Naim Banchereau Dolinay Parnell Dolinay Parnell Lill Liesenfeld Rogers Liesenfeld Giamarellos-Bourboulis	Sepas. mary cases from burkholderia Pendenkon and the separation of the separation pendenkon separation of the separation Pendemic HINI in CU Bacterial Infection or other inflammatory conditions Secret HINI Separation and the separation of the separation Secret HINI Separation and the separation of the separation for the second second second second second reference second second second second second reference second second second second second Patients undergoing mechanical ventilation with second second second second second second pendential second second second second second second pendential second second second second second second second pendential second second second second second second pendential second second second second second second pendential second second second second second second second pendential second s	55 20 14 274 14 21 11 80 81 44 54 123 52 21 59 13 3	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 56.0 (1QR 42.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 57.0 (1QR 29.0-56.0) 57.0 (1QR 24.0-72.0) 65.0 (1QR 37.0-85.0) 72.0 (1QR 71.0-57.40.0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62) 2 (67)	IIIumina IIIumina IIIumina Affymetrix Affymetrix Affymetrix IIIumina IIIumina IIIumina IIIumina IIIumina Affymetrix NanoString® NanoString NanoString NanoString	Thailand Thailand Unk. Nerway Australia Unk. Unk. United States of America United States of America Australia Vietnam United States of America Estonia United States of America United States of America USA/Greece USA	$\begin{array}{c} 45(82)\\ 15(75)\\ 0(0)\\ 52(19)\\ 0(0)\\ 10(48)\\ 0(0)\\ 0(0)\\ 59(73)\\ 0(0)\\ 16(30)\\ 0(0)\\ 8(15)\\ 21(100)\\ 24(41)\\ 6(46)\\ 5(38)\\ 3(100)\\ \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 0 (0) 0 (0) 0 (0) 22 (42) 0 (0) 7 (12) 1 (8) 8 (62) 0 (0)	$\begin{array}{c} 10 \ (18) \\ \hline 10 \ (18) \\ \hline 5 \ (25) \\ \hline 4 \ (28) \\ \hline 222 \ (81) \\ \hline 7 \ (50) \\ \hline 11 \ (52) \\ \hline 0 \ (0) \\ \hline 22 \ (27) \\ \hline 44 \ (100) \\ \hline 30 \ (56) \\ \hline 0 \ (0) \\ \hline 22 \ (42) \\ \hline 0 \ (0) \\ \hline 28 \ (47) \\ \hline 6 \ (46) \\ \hline 0 \ (0) \\ \hline 0 \ (0) \\ \hline \end{array}$
CSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE21802 GSE22998 GSE29385 GSE29385 GSE29385 GSE29385 GSE29391 GSE29395 GSE29395 GSE29395 GSE29395 GSE29395 GSE301996 GSE40035 INF-02 INF-02 INF-02 INF-015-01 INF-03-04	Pankla Penkla Bernojo-Martin Bery Berdal Sutherland Naim Banchereau Dolinay Parnell Nguyen Hu Lall Liesenfeld Rogers Liesenfeld Giamarellos-Bourboulis Ieleviss	Sepas, mary cases from burkholderia Sepais, mary cases from Pendemie HINI in ICU Bacieral infection or other inflammatory conditions Severe HINI Sepaie and Severe HINI Severe HINI Patients with output Bacterial or influenza Patients undergoing mechanical ventilation History and SIRS Febril e children with decrease Febrile Severe Sinter Severe Severe History tract infections ED and hospitalized sepais Patients with Bacterial upper respiratory tract infections ED and hospitalized sepais Separet hock registratory to ICI	55 20 14 274 121 11 80 81 44 54 123 52 21 59 13 3 6	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 6.5 (1QR 1.92-11.0) 56.0 (1QR 45.0-70.0) 12.0 (1QR 10.0-14.0) 57.0 (1QR 45.0-70.5) 39.0 (1QR 45.0-70.5) 39.0 (1QR 29.0-56.0) 57.0 (1QR 73.0-74.0) 65.0 (1QR 73.0-85.0) 72.0 (1QR 70.5-74.0) 64.5 (1QR 60.05.71 0)	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62) 2 (67) 3 (50)	Illumina Illumina Illumina Afymetrix Afymetrix Illumina Illumina Illumina Illumina Illumina Illumina Aflymetrix NanoString NanoString NanoString	Thailand Thailand Unk. Unk. Unk Unk Unk United States of America Australia Vietnam United States of America Estonia UsA/Greece USA Greece IISA	$\begin{array}{c} 45 (82) \\ 15 (75) \\ 0 (0) \\ 52 (19) \\ 0 (0) \\ 10 (48) \\ 0 (0) \\ 0 (0) \\ 59 (73) \\ 0 (0) \\ 59 (73) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 16 (30) \\ 0 (0) \\ 5 (30) \\ 16 (30) \\ 10 (30) \\ 16 (30) \\ 10 (30) \\ $	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 8 (15) 123 (100) 0 (0) 8 (15) 123 (100) 7 (12) 1 (8) 8 (62) 0 (0) 0 (0)	$\begin{array}{c} 10(18) \\ 5(25) \\ 4(28) \\ 7(50) \\ 7(50) \\ 0(0) \\ 0(0) \\ 22(27) \\ 44(100) \\ 30(55) \\ 0(0) \\ 22(47) \\ 6(46) \\ 0(0) \\ 28(47) \\ 6(46) \\ 0(0$
GSE13015gpl6102 GSE13015gpl6047 GSE21802 GSE21802 GSE22098 GSE22098 GSE2985 GSE2985 GSE30119 GSE2985 GSE30119 GSE30119 GSE4015 GSE40165 GSE40165 GSE40165 GSE40165 GSE40180 INF-02 INF-118-01 INF-118-03 INF-118-04 INF-118-01	Pankla Pankla Bernejo-Martin Berty Bertal Sutherland Naim Naim Dolinay Panell Dolinay Panell Nguyen Hu Lill Lill Liesenfeld Rögers Liesenfeld Giamarellos-Bourboulis Iglesias	Sepas. mary cases from burkhoder: Pandemic HIN in CU Bacterial Infection or other inflammatory conditions Severe HINI Severe H	55 20 14 274 14 21 11 80 81 44 52 21 59 13 3 6	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 13.5-54.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 12.0-20.0) 50.0 (1QR 42.0-59.0) 50.0 (1QR 45.0-59.0) 50.0 (1QR 45.0-59.0) 50.0 (1QR 43.0-79.0) 50.0 (1QR 13.0-14.0) 0.92 (1QR 0.3.3-1.6) 57.0 (1QR 33.0-70.5) 57.0 (1QR 42.0-72.0) 65.0 (1QR 37.0-85.0) 72.0 (1QR 70.5-74.0) 64.5 (1QR 60.25-71.0)	$\begin{array}{c} 26 \ (47) \\ 6 \ (30) \\ 0 \ (0) \\ 59 \ (22) \\ 10 \ (71) \\ 0 \ (0) \\ 48 \ (60) \\ 34 \ (42) \\ 13 \ (30) \\ 20 \ (37) \\ 85 \ (69) \\ 17 \ (33) \\ 0 \ (0) \\ 17 \ (28) \\ 10 \ (77) \\ 8 \ (62) \\ 2 \ (67) \\ 3 \ (50) \end{array}$	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina Illumina NanoString [®] NanoString [®] NanoString NanoString NanoString NanoString	Thailand Thailand Unk. Norway Australia Unk. Unk. United States of America States of America Australia Vietnam United States of America Estonia USA/Greece USA	$\begin{array}{c} 45 \ (82) \\ 15 \ (75) \\ 0 \ (0) \\ 52 \ (19) \\ 0 \ (0) \\ 10 \ (48) \\ 0 \ (0) \\ \hline \end{array} \\ \begin{array}{c} 60 \ (0) \\ 59 \ (73) \\ 0 \ (0) \\ 16 \ (30) \\ 0 \ (0) \\ 16 \ (30) \\ 0 \ (0) \\ 24 \ (41) \\ 6 \ (46) \\ 5 \ (38) \\ 3 \ (100) \\ 6 \ (100) \\ \end{array}$	0 (0) 0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 80 (100) 0 (0) 0 (0) 20 (10) 0 (0) 20 (0) 123 (100) 22 (42) 0 (0) 7 (12) 1 (8) 8 (62) 0 (0) 0 (0)	$\begin{array}{c} 10(18) \\ 5(25) \\ 4(28) \\ 7(50) \\ 0(0) \\ 0(0) \\ 22(27) \\ 44(100) \\ 30(56) \\ 0(0) \\ 22(42) \\ 0(0) \\ 22(42) \\ 0(0) \\ 6(46) \\ 0(0) \\ 0(0) \\ \end{array}$
CSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22098 GSE27331 GSE29385 GSE29385 GSE29385 GSE29391 GSE29385 GSE23707 GSE3019 GSE3017 GSE40196 GSE40196 GSE40396 INF-02 INF-02 INF-01 INF-01 INF-03 INF-03 INF-04 IN	Pankla Penkla Bernojo-Martin Berny Berdal Sutherland Naim Banchereau Dolinay Pamell Nguyen Ha Lill Liesenfeld Rogers Liesenfeld Giamarellos-Bourboulis Bauer	Sepas, mary cases from burkholderia Sepais, mary cases from Pendemic HIN in a CU Bacieral Infection or other inflammatory conditions Severe HINI Sepais or post-surgical SIRS Patients with surgected dengue Febrile patients with oral infection, mostly influenza Pediatric patients with oral patients with community-acquired Bacterial or influenza A prediation sepais or SIRS on element Bacterial or influenza A presentation or SIRS Petrile children with dengue Febrile children with dengue Febrile children with dengue Febrile children with dengue Febrile children with bacterial community acquired or single sepais of SIRS Petrile children with bacterial common fish factor for acute respiratory distress syndrome Patients with bacterial upper respiratory tract infections ED and hospitalized sepais Dath bacterial gepsis patients with bacterial gepsis patients with bacterial gepsis patients with activities of the sepais of septic shock patients with activities of the septies of septic shock patients with activities of the septies of septies bacterial Septis of septic shock patients admitted to ICU Patients with activities of the septies of septies hock patients with activities of the septies of septies hock patients admitted to ICU	55 20 14 21 11 80 81 44 52 21 52 21 59 13 3 6 25	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9-21.10.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-70.0) 12.0 (1QR 10.0-14.0) 0.92 (1QR 03.31.6) 57.0 (1QR 73.0-70.5) 39.0 (1QR 29.0-56.0) 57.0 (1QR 73.0-85.0) 72.0 (1QR 65.0-87.0) 64.5 (1QR 60.2-71.0) 64.5 (1QR 60.2-71.0) 65.0 (1QR 70.5-74.0) 64.5 (1QR 60.2-71.0) 65.0 (1QR 60.2-71.0) 64.5 (1QR 60.2-71.0) 65.0 (1QR 60.2-71.0) 65.0 (1QR 60.2-71.0) 64.5 (1QR 60.2-71.0) 65.0 (1QR 60.2-71.0) 64.5 (1QR 60.2-71.0) 65.0 (1QR 60.2-71.0) 65.0 (1QR 60.2-71.0) 64.5 (1QR 60.2-71.0) 65.0 (1QR 70.2-71.0) 65.0 (1QR	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62) 2 (67) 3 (50) 14 (55)	Illumina Illumina Illumina Affymetrix Affymetrix Illumina Illumina Illumina Illumina Illumina Illumina Affymetrix NanoString NanoString NanoString	Thailand Thailand Unk Norway Australia Unk Unk Unk Unk Unk States of America Australia Vietnam United States of America Estonia USA/Greece USA Greece USA	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0) 8 (15) 21 (100) 8 (15) 21 (100) 6 (46) 5 (38) 3 (100) 6 (100)	$\begin{array}{c} 6 \ (6) \\ \hline 0 \ (0) \\ 10 \ (71) \\ \hline 0 \ (0) \\ 7 \ (50) \\ \hline 0 \ (0) \\ \hline 11 \ (100) \\ 80 \ (100) \\ \hline 0 \ (0) \\ \hline 0 \ (0) \\ \hline 0 \ (0) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 23 \ (100) \\ \hline 123 \ (100) \\ \hline 11 \ (100) \ (100) \\ \hline 11 \ (100) \ (100) \ (100) \\ \hline 11 \ (100) \ ($	$\begin{array}{c} 10(18) \\ 5(25) \\ 4(28) \\ 7(50) \\ 0(0) \\ 222(27) \\ 11(52) \\ 0(0) \\ 22(27) \\ 44(100) \\ 30(56) \\ 0(0) \\ 22(47) \\ 6(46) \\ 0(0) \\ 28(47) \\ 6(46) \\ 0(0) \\ $
CSE13015gpl6102 GSE13015gpl6947 GSE21802 GSE22998 GSE27181 GSE2985 GSE30119 GSE3950 GSE30119 GSE32707 GSE40012 GSE400000	Pankla Bernejo-Martin Berry Berdal Sutherland Naim Naim Dolinay Pamell Nguyen Hu Linesenfeld Liesenfeld Giamarello-Bourboulis Iglesias Bauer	Sepas, mary cases from burkhoder, Pendemie HINI in CU Pandemie HINI in CU Bacieral infection or other inflammatory conditions Severe HINI Seyes or post-surgical SIRS Patients with suspected dengue referite panets with community-acquired bacterial infection Patients undergoing mechanical ventilation with seysis or SIRS or neither Bacterial or influenza A patients with community-acquired bacterial infection Patients undergoing mechanical ventilation with seysis or SIRS or neither Bacterial or influenza A presention of SIRS as reference on SIRS as reference on SIRS as reference on the second or viral infections with saccinal or viral infections with Patients with suspected infections/ Patients with sub-reference on Patients with bacterial appear reprintery text infections ED and hospitalized sepsis patients Sepsis or septies shock patients with suppected infections/ patients with suppected infections/ patients with suppected infections/ patients with suppected infections/ Sepsis or septies shock patients with suppected infections/ Sepsis or septies shock patients with suppected infections/ patients with suppected infections/ Septis or septies thock patients with suppected infections/ Septis or septies thock patients with suppected infections/ Septis or septies and compare Patients with suppected infections/ Septis or septies and compare patients and septient of course	55 20 14 274 14 21 11 80 81 44 54 123 52 21 59 13 13 3 6 25	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.92-11.0) 50.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-50.0) 57.0 (1QR 37.0-85.0) 57.0 (1QR 37.0-85.0) 57.0 (1QR 37.0-85.0) 57.0 (1QR 37.0-85.0) 57.0 (1QR 40.2-77.0) 65.0 (1QR 40.2-77.0) 60.0 (1QR	26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62) 2 (67) 3 (50) 14 (56)	Illumina Illumina Illumina Allymetrik Allymetrik Illumina Illumina Illumina Illumina Illumina Illumina Illumina Illumina Illumina Nanostring NanoString NanoString NanoString	Thailand Thailand Unk. Nervey Australia Unk. Unk. United States of America Australia Vietnam United States of America Australia Vietnam United States of America Usa Usa Greece USA Greece USA	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30) 0 (0) 16 (30) 0 (0) 24 (41) 6 (46) 5 (38) 3 (100) 16 (64)	0 (0) 10 (71) 0 (0) 7 (50) 0 (0) 11 (100) 80 (100) 0 (0) 0 (0) 8 (15) 123 (100) 22 (42) 0 (0) 22 (42) 0 (0) 7 (12) 1 (8) 8 (62) 0 (0) 5 (20)	$\begin{array}{c} 10(18) \\ 5(25) \\ 4(28) \\ 7(50) \\ 11(52) \\ 10(52) \\ 11(52) \\ 10(52) \\ 11(52) \\ 10(52) \\ 11(52) \\ 10$
CSE13015gp16102 GSE13015gp16947 GSE21802 GSE22098 GSE22098 GSE22998 GSE29981 GSE29991 GSE29991 GSE29991 GSE29991 GSE20707 GSE40012 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE40002 GSE4002 GSE4002 GSE4002 GSE4002 GSE4002 GSE4002 GSE4002	Pankla Pankla Berncjo-Martin Berdal Sutherland Naim Banchereau Dolinay Pamell Nguyen Hu Lill Liesenfeld Rogers Liesenfeld Giamarellox-Bourboulis Iglesias Bauer	Sepas, may cases from burkholderia Sepais, may cases from burkholderia Bacieral Infection or other infammatory conditions Severe HINI Sepais or post-augical SIRS Patients with suspected dengue Febrile patients with rommunity-acquired bacterial infection, mostly influenza Pediatre patients with community-acquired bacterial infection Patients unit community-acquired bacterial infection Patients with community-acquired bacterial infection with regroupin grBS or earlier Bacterial or influenza A precumoia or SIRS Pedrile children with dengue Pedrile statistics with saterial community-acquired bacterial meningitis Patients with saterial autorespis Patients with saterial upper respiratory tract infections ED and hospitalized sepsis patients Sepsis or septis shock patients with saypected infection/spiss Patients with saypected infection/spiss	55 20 14 274 14 21 11 80 81 44 52 21 59 13 3 6 25 5	54.0 (1QR 48.0-61.0) 49.0 (1QR 43.5-59.5) Unk. 16.0 (1QR 11.0-26.0) 39.0 (1QR 35.75-44.75) Unk. 25.0 (1QR 22.0-40.0) 6.5 (1QR 1.9-21.10.0) 56.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-59.0) 59.0 (1QR 45.0-57.0) 12.0 (1QR 10.0-14.0) 0.92 (1QR 33.0-70.5) 39.0 (1QR 45.0-70.5) 57.0 (1QR 70.5-74.0) 64.5 (1QR 60.25-71.0) 77.0 (1QR 65.0-86.0) 51.0 (1QR 70.5-74.0) 51.0 (26 (47) 6 (30) 0 (0) 59 (22) 10 (71) 0 (0) 48 (60) 34 (42) 13 (30) 20 (37) 85 (69) 17 (33) 0 (0) 17 (28) 10 (77) 8 (62) 2 (67) 3 (50) 14 (56) 5 (100)	Illumina Illumina Illumina Aflymetrix Aflymetrix Illumina Illumina Illumina Illumina Illumina Illumina Aflymetrix NanoString NanoString NanoString NanoString	Thailand Thailand Unk Unk Unk Norway Australia Unk Unk United States of America USA/Greece USA Greece USA	45 (82) 15 (75) 0 (0) 52 (19) 0 (0) 10 (48) 0 (0) 59 (73) 0 (0) 16 (30), 0 (0) 8 (15) 21 (100) 8 (15) 21 (100) 8 (46) 5 (38) 3 (100) 6 (100) 16 (64)	$\begin{array}{c} 6 \ (6) \\ \hline 0 \ (0) \\ 10 \ (71) \\ \hline 0 \ (0) \\ 7 \ (50) \\ \hline 0 \ (0) \\ \hline 11 \ (100) \\ 80 \ (100) \\ \hline 0 \ (0) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 123 \ (100) \\ \hline 22 \ (42) \\ \hline 11 \ (8) \\ \hline 8 \ (62) \\ \hline 0 \ (0) \\ \hline 5 \ (20) \\ \hline 5 \ (20) \\ \hline 5 \ (20) \\ \hline \end{array}$	$\begin{array}{c} 10(18) \\ 5(25) \\ 4(28) \\ 222(81) \\ 7(50) \\ 0(0) \\ 0(0) \\ 22(27) \\ 44(100) \\ 30(56) \\ 0(0) \\ 22(42) \\ 0(0) \\ 28(47) \\ 6(46) \\ 0(0) \\ 28(47) \\ 6(46) \\ 0(0) \\ $

Study	Author	Description	N	A ga*	Mala ⁺	Platform	Country	30 day Mortality	Survival
Study	Autioi	Patiants hospitalized with COPD	18	Age	wate	Flatiofin	Country	50-day Mortanty	Survivar
EMEVD2590	Almonos	r attents hospitalized with COLD	27	Unit	16 (50)	Anilant	Sauia	0.00	27 (100)
EMEAF3309	Almansa	Surgical patients with capsic	21	Ulik.	10 (39)	Agnem	Span	0(0)	27(100)
EMTAD1549	Almonos	(EVDEESS)	140	72.0 (IOP 61.0 78.0)	05 (69)	Anilant	Sauia	12(0)	128 (01)
EMIADIJ46	Aimailsa	(EAFRESS)	140	72.0 (IQR 01.0-78.0)	95 (08)	Agnent	Span	12 (9)	128 (91)
EMTAB5102	van den Ham	Patients with dengue	21	20.0 (IQR 17.0-28.0)	11 (52)	Anymetrix	Indonesia	0 (0)	21 (100)
CSE42026	Harbarg	bacterial infaction	50	1.25 (TOP 0.38 4 0)	33 (56)	Illumina	United Kingdom	0.00	59 (100)
CSE47655	Stone	Savara anaphylavis in ED	6	1.2.5 (IQR 0.56-4.0)	0.00	Affymatrix	Australia	0(0)	6(100)
CSE51909	Kasisas	Detients with densus	20	Unk.	0(0)	Affyniculx	Theiland	0(0)	28 (100)
GSEJ1808	Kwissa	Fatients with deligue	20	Ulik.	0(0)	Anymeurx	Thanand	0(0)	28 (100)
GSE54514	Parneli	Sepsis patients in ICU	19	Unk.	0(0)	Illumina	Unk.	4 (21)	15 (79)
GSE37183	Shellor	Februe cunidien with SJIA	11	3.0 (IQK 3.3-7.3)	3 (43)	Inumna	United States of America	0(0)	11(100)
GSE25504	Smith	Neonatal sepsis	33	0.0 (IQR 0.0-0.0)	19 (57)	Anymetrix	United Kingdom	2 (0)	31 (94)
GSE00244	Suarez	Lower respiratory tract infections	118	63.0 (IQR 50.0-77.0)	37 (31)	Illumina	United States of America	0 (0)	118 (100)
00000000		retiatric patients with Kawasaki		A 85 (COD 4 10 4 05)	100.000			0.00	
GSE63881	Hoang	disease	1/1	2.75 (IQR 1.42-4.25)	102 (60)	Illumina	United States of America	0 (0)	171 (100)
CEEGAASG	Mahaian	reprise infants with bacterial	200	0.1 (TOP 0.06 0.12)	106 (52)	Thuming	United States of America	0.00	200 (100)
03E04430	wanajan	or vital infection in ED	200	0.1 (IQK 0.00-0.13)	100 (33)	mumma	United States of America	0(0)	200 (100)
CODUCTO	6.1.1	Suspected but negative for	100	11-1	0.00	100	11-1	22 (22)	02 (70)
G5E05082	Sciciuna	CAP Padiatric [C1] patiante with	106	Unk.	0(0)	Anymetrix	Unk.	23 (22)	85(78)
00000000	6	remaine ico patients with	220	2.45 (COD 1.0.5.00)	04.7410	100	The LOW CARD	20 (12)	201 (00)
GSE00099	Sweeney	sepsis, SIRS	229	2.45 (IQR 1.0-5.88)	94 (41)	Anymetrix	United States of America	28 (12)	201 (88)
G5E07059	Heinonen	Children with rninovirus	80	0.85 (IQK 0.5-1.29)	33 (00)	iiiumina	United States of America	0(0)	80 (100)
		Outpatients with acute respiratory							
GSE68310	Zhai	viral infections	104	20.96 (IQR 20.09-22.76)	50 (48)	Illumina	United States of America	0 (0)	104 (100)
		Pediatric patients with bacterial							
GSE72810	Herberg	or viral infection	72	1.83 (IQR 0.88-3.29)	8(11)	Illumina	United Kingdom	0 (0)	72 (100)
		Pediatric patients with bacterial							
GSE73461	Wright	or viral infection or other inflammatory diseases	404	2.79 (IQR 0.92-8.81)	165 (41)	Illumina	United Kingdom	0 (0)	404 (100)
		Infants with respiratory synctial							
GSE77087	de Steenhuijsen Piters	virus	41	0.45 (IQR 0.14-0.69)	25 (61)	Illumina	United States of America	0 (0)	41 (100)
00000000		Moderate and severe influenza			10 10 1			0.00	
GSE82050	lang	intection	24	64.5 (IQR 48.5-74.25)	13 (54)	Agilent	Germany	0 (0)	24 (100)
GSE95255	venet	Septic snock	51	00.0 (IQR 55.5-75.5)	31 (01)	Anymetrix	France	17 (55)	34 (07)
00000000		Infants with respiratory synctial	10	0.05 (COD 0.18 0.10)	20. (62)			0.00	1000
GSE103842	Rodriguez-Fernandez	virus	62	0.25 (IQR 0.17-0.44)	39 (63)	Illumina	United States of America	0 (0)	62 (100)
		Sepsis, many cases from		54.0 COD 10.0 CL 0				15.00	10.000
GSE13015gpl6102	Pankla	burkholderia	55	54.0 (IQR 48.0-61.0)	26 (47)	Illumina	Thailand	15 (27)	40 (73)
		Sepsis, many cases from		10.0 200 10.5 50 5	6.000				0.000
GSE13015gpl6947	Pankla	burkholderia	15	49.0 (IQR 43.5-59.5)	6 (40)	Illumina	Thailand	7 (47)	8 (53)
COTRADOOR	D	Bacterial infection or other	102	160.000 110.000	60 (21)	manter	11-1	0.00	102 (100)
GSE22098	Berry	innammatory conditions	195	16.0 (IQR 11.0-26.0)	59 (31)	Illumina	Unk.	0(0)	195 (100)
GSE2/151	Berdai	Severe HINI	/	38.0 (IQR 33.0-50.0)	0 (80)	Anymetrix	INOFWAY	2 (28)	5(71)
GSE28991	INaim	Patients with suspected dengue	11	Unk.	0(0)	Illumina	Unk.	0 (0)	11 (100)
COTTONIA	Destaura	reutatric patients with community-acquired	50	6.6.000 1.00 11.00	24 (877)	manta	The low of the state of the sta	0.00	E0 (100)
GSE30119	Banchereau	bacterial infection	59	0.5 (IQK 1.92-11.0)	54 (57)	Illumina	United States of America	0 (0)	59 (100)
00000000	DUT	rations undergoing mechanical ventilation	60		12 (10)	manta	The low of the state of the sta	25.000	11/60
GSE32707	Dolinay	with sepsis or SIRS or neither	69	20.0 (IQK 45.0-59.0)	13 (19)	Illumina	United States of America	25 (36)	44 (64)
CSE40012	Parnall	Bacteriai or influenza A	20	59.0 (TOP 46.5.67.0)	20 (51)	Illumine	Australia	5 (12)	34 (97)
GSE40012	Nauvan	Fabrila children with dengue	122	120 (IQK 40.3-07.0)	20 (31)	Illumina	Vietnam	0(0)	123 (100)
03040103	nguyen	Febrile children with bacterial	123	12.0 (IQK 10.0-14.0)	65 (09)	mumma	victiani	0(0)	123 (100)
GSE40396	Hu	or viral infection	30	0.92 (TOR 0.33-1.6)	17 (56)	Illumina	United States of America	0.00	30 (100)
GSE40586	Lill	Community-acquired bacterial meningitis	21	57.0 (IOR 53.0-70.5)	0.00	Affymetrix	Estonia	2(10)	19 (90)
001040000		Patients with suspected infection/sensis	~1	57.0 (121 55.0-70.5)	0 (0)	7 my meurx	Latolila	2 (10)	
INE-02	Liesenfeld	in ED	61	41.0 (IOR 29.0-57.0)	18 (30)	NanoString®	USA/Greece	2 (3)	59 (97)
111 02	Elesenteid	ICU patients with at	01	41.0 (1012).0 57.0)	10 (50)	runooung	USIVGICCC	2 (3)	37(71)
INF-IIS-01	Rogers	least one risk factor for acute respiratory distance condrome	11	57.0 (IOR 39.5-81.5)	10 (91)	NanoString	USA	5 (45)	6 (55)
101 113-01	nogers	Patients with bacterial upper		57.0 (121 59.5-61.5)	10 (91)	runosung	C3A	5 (45)	0 (33)
INF-IIS-02	Liesenfeld	respiratory tract infections	19	77.0 (IOR 51 5-85.0)	12 (63)	NanoString	Greece	9 (47)	10 (53)
141-113-03	Liesenieu	ED and hospitalized sensis	1.7		12 (03)	. ano.5u mg	Gitte	2(77)	10(55)
INF US 04	Giamarallos Rourboulis	patients	2	72.0 (TOP 70.5.74.0)	2 (67)	NanoString	Graaca	1 (33)	2 (67)
1111-110-04	communerer of the second of th	Sensis or sentic shock	5	72.0 (IQIC 10.3-14.0)	2 (07)	. ano.5u mg	Gitte	1 (55)	2(07)
INFE US 10	Interior	notiants admitted to ICU	6	64.5 (TOP 60.25.71.0)	3 (50)	NanoString	USA	1 (17)	5 (82)
101-113-10	igiesias	Patients admitted to reco	U	04.5 (IQK 00.25-71.0)	5 (50)	ranosung	USA	1(17)	5 (65)
INF US 11	Ramar	in ED	25	77.0 (10) 65.0.86.0	14 (56)	NanoString	Garmany	4(16)	21 (84)
1007-110-11	Dauci	III ED Patiante with suspected severe	2.5	77.0 (IQK 03.0-80.0)	14 (30)	ranosung	Germany	+(10)	21(04)
INTE THE OF	A 1	rations with suspected severe	¢	810 (OP 770 82 0)	5 (100)	None State	Sauia	0.00	5 (100)
nvr-115-21	Annansa	respiratory infection	5	01.0 (IQK //.0-82.0)	5 (100)	ivanoString	opain	0(0)	5 (100)

Table 4: **Development datasets for prognostic classifier.** SJIA: Systemic Juvenile Idiopathic Arthritis.

8.3 Feature selection

Feature selection is illustrated in Fig. 5. It was assay-guided in the sense that all selected features must have exhibited acceptable amplification by LAMP.

- 1. The "benchtop LAMP" experiments revealed lack of amplification of 10 genes, which were replaced by 10 genes with a good trade-off between the LAMP amplification and clinical performance [23].
- 2. Subsequent experiments on the target platform revealed lack of amplification of additional 11 genes, due to differences between the benchtop and Myrna[™] amplification. The non-amplified genes were replaced with genes which had good trade-off between Myrna amplification and clinical performance using the process shown in Fig. 5. Briefly, we started with a new sorted list of predictive and prognostic genes identified using midplex-genomewide search for biological signal. We then went down the list and synthesized Myrna primers (chemicals used to detect mRNA), where possible. The highest ranked gene with amplified primer replaced the highest ranked gene without the amplification in the existing signature. The ranking was based on a heuristic function comprising Shapley value for the gene and concordance with the NanoString platform. The procedure was repeated until all non-amplified genes were replaced.
- 3. Following the transfer of the gene set to Myrna, we attempted reduction of the signature size using a greedy forward search.

The final set of genes chosen by the feature selection process is shown in Fig. 6.



Figure 5: **Overview of selection of gene features compatible with LAMP amplification.** GWS: Genome-wide search; MA: microarray; NS: NanoString[®].

Table 5: **Impact of CF on BVN classifier accuracy and concordance.** The CF experiments used the filtering value of 0.92 Pearson correlation coefficient between probabilities assigned to replicate samples.

Algorithm	Statistic	With CF	Without CF
LOGR	CV mAUROC	0.886	0.885
LOGR	Validation mAUROC	0.834	0.833
LOGR	Bacterial PCC	0.938	0.939
LOGR	Viral PCC	0.962	0.963
MLP	CV mAUROC	0.895	0.897
MLP	Validation mAUROC	0.832	0.834
MLP	Bacterial PCC	0.942	0.891
MLP	Viral PCC	0.960	0.940

8.4 Concordance Filtering ablation

To estimate the impact of CF, we performed an ablation study by evaluating performance of the leading classifiers (MLP and LOGR) with and without the CF. Table 5 shows that CF significantly improves the PCC for MLP classifier, without loss of accuracy.



Figure 6: **TriVerity[™] genes, ordered by absolute values of Shapley values for each probability. a,** BVN classifier. **b,** SEV classifier.