M4CXR: Exploring Multi-task Potentials of Multi-modal Large Language Models for Chest X-ray Interpretation

Jonggwon Park, Soobum Kim, Byungmu Yoon, Jihun Hyun, Kyoyun Choi*

DEEPNOID Inc., Seoul, Korea {jgpark, soobumk, bmyoon, jhyun, kychoi}@deepnoid.com

Abstract

The rapid evolution of artificial intelligence, especially in large language models (LLMs), has significantly impacted various domains, including healthcare. In chest X-ray (CXR) analysis, previous studies have employed LLMs, but with limitations: either underutilizing the multi-tasking capabilities of LLMs or lacking clinical accuracy. This paper presents M4CXR, a multi-modal LLM designed to enhance CXR interpretation. The model is trained on a visual instructionfollowing dataset that integrates various task-specific datasets in a conversational format. As a result, the model supports multiple tasks such as medical report generation (MRG), visual grounding, and visual question answering (VQA). M4CXR achieves state-of-the-art clinical accuracy in MRG by employing a chain-of-thought prompting strategy, in which it identifies findings in CXR images and subsequently generates corresponding reports. The model is adaptable to various MRG scenarios depending on the available inputs, such as single-image, multi-image, and multi-study contexts. In addition to MRG, M4CXR performs visual grounding at a level comparable to specialized models and also demonstrates outstanding performance in VQA. Both quantitative and qualitative assessments reveal M4CXR's versatility in MRG, visual grounding, and VQA, while consistently maintaining clinical accuracy.

Introduction

Recent advancements in artificial intelligence, particularly large language models (LLMs), have led to their widespread application across various fields, including healthcare. Numerous studies are exploring diverse methods of application for improving healthcare outcomes, such as personalized treatment plans, clinical decision support systems, and enhancing medical education (Zheng et al. 2024).

Within the medical domain, this paper focuses on chest X-ray (CXR) interpretation. Although many studies have investigated LLMs for CXR, most encounter one of two issues: they either underutilize the capabilities of LLMs or struggle with ensuring clinical accuracy. Despite LLMs' ability to perform various tasks through conversational interaction, many studies focus on a single task, typically medical report generation (MRG) (Zhou et al. 2024; Hyland et al. 2024; Chaves et al. 2024). While some studies effectively

leverage LLMs to comprehend questions and generate appropriate responses (Thawkar et al. 2023; Pellegrini et al. 2023; Lee et al. 2024; Chen et al. 2024), the clinical accuracy of these conversational outcomes, including generated reports, falls short of expectations.

In this work, we propose M4CXR, an LLM for CXR adept at handling four "Multi" aspects: Multi-modal, Multitask, Multi-image input, and Multi-turn chatting. We train the model on a dataset constructed by integrating various task-specific datasets, enabling M4CXR to excel in MRG, image understanding, and visual question answering (VOA) tasks, as depicted in Figure 1. In MRG, M4CXR improves clinical accuracy by leveraging LLMs' reasoning abilities through a chain-of-thought (CoT) (Wei et al. 2022) reasoning, where the model identifies findings in the image and generates reports based on these results. Additionally, we conduct experiments across various scenarios with inputs extending beyond a single image, including CXR images from different views and prior patient studies. Besides MRG, M4CXR demonstrates strong multi-tasking capability with successful application in visual grounding and VQA, showcasing its adaptability in diverse clinical contexts.

In summary, our main contributions are as follows:

- We propose M4CXR, an MLLM designed for CXR interpretation, capable of handling multiple tasks. To enable multi-tasking, we assemble a visual instructionfollowing dataset from diverse CXR tasks.
- By adopting a novel CoT reasoning process, M4CXR achieves state-of-the-art clinical accuracy in CXR report generation.
- M4CXR effectively utilizes multiple images and reports, allowing for its applicability across different scenarios.
- Beyond MRG, M4CXR demonstrates remarkable performance in visual grounding and VQA.

Related Works

Multi-modal Large Language Models

As research on LLMs progresses, their advantages continue to emerge. It is widely known that LLMs are capable of multi-task learning (Raffel et al. 2020), and Wei et al. (2022) showed that LLMs can reason through multi-turn chats using

^{*}Corresponding author



Figure 1: Overview of the multi-tasking capabilities of M4CXR. Facilitated by CoT prompting in MRG, M4CXR produces clinically accurate reports and adapts to various scenarios. Additionally, M4CXR can ground the locations described in the report or answer questions based on chest X-ray images.

CoT prompting. The rise of multi-modal LLMs, which utilize the long context length of LLMs to process visual content, further amplifies these advantages. MLLMs can handle tasks that require comprehension of spatial regions, and visual instruction tuning (Liu et al. 2024) allows them to answer free-form questions about images.

For MLLMs, bridging the embedding spaces of vision and language is crucial. A common approach is to freeze the pretrained vision model and LLM, and train only the bridging module (Alayrac et al. 2022; Zhu et al. 2023), also known as the projector. The projector's structure varies, from a linear layer to Q-Former (Li et al. 2023a), C-Abstractor (Cha et al. 2024), and more. Drawing insights from these works, we employ a two-phase training strategy: pre-training the projector first, and then fine-tuning the entire model with visual instruction tuning.

Chest X-ray Report Generation

CXR report generation via deep learning has been extensively studied over time (Shin et al. 2016; Jing, Xie, and Xing 2018; Zhang et al. 2020; Tanida et al. 2023). Since the advent of LLMs, training strategies such as freezing pretrained models (Thawkar et al. 2023) or visual instruction tuning (Lee et al. 2024) have been equivalently implemented in the CXR domain. However, a significant drawback of these models is the lack of clinical accuracy in the generated reports. Efforts to address these shortcomings include adopting a separate disease classification module (Jin et al. 2024; Zhou et al. 2024), increasing the visual encoder's resolution (Chaves et al. 2024), or harnessing additional input (Bannur et al. 2024). In this work, we take a novel yet simple approach: leveraging the reasoning abilities of LLMs.

Multi-Tasking in Chest X-ray Interpretation

Confining the use of LLMs solely to report generation underutilizes their potential. MLLMs that capitalize on the characteristics of LLMs to enable conversations based on CXR images include XrayGPT (Thawkar et al. 2023), Ra-Dialog (Pellegrini et al. 2023), and LLM-CXR (Lee et al. 2024). CheXagent (Chen et al. 2024) is a multi-tasking CXR foundation model trained on various tasks, similar to our approach. Yet, the clinical accuracy of these models' responses remains suboptimal.

Medical foundation models such as Med-Gemini (Yang et al. 2024) and MedPaLM-M (Tu et al. 2024) support CXR interpretation with conversational capabilities, showing satisfactory clinical accuracy. However, these models are not trained on spatial cognition tasks like visual grounding, which involves identifying regions in an image corresponding to a phrase in the text. For CXR, Boecking et al. (2022a) released the benchmark dataset MS-CXR for phrase grounding, and MedRPG (Chen et al. 2023b) specializes in identifying bounding boxes for phrases in CXR images. The latest concurrent work by Bannur et al. (2024) introduced MAIRA-2, an MRG model for CXR capable of grounding but limited to grounded report generation. In contrast, M4CXR handles multiple tasks such as MRG, VQA, and visual grounding, while maintaining clinical accuracy.

Methods

M4CXR

Architecture. Figure 2 (a) illustrates the overall architecture of M4CXR. Following LLaVA (Liu et al. 2024), M4CXR includes a vision encoder, a projector, and an LLM, denoted as $\mathcal{M}_V, \mathcal{M}_P$, and \mathcal{M}_L , respectively. Given n CXR images x_1, \ldots, x_n , the vision encoder converts each image x_i into visual feature f_i^v , which the projector then transforms into a sequence of visual token embeddings X_i^v :

$$X_i^v = \mathcal{M}_P(f_i^v) = \mathcal{M}_P(\mathcal{M}_V(x_i)), \quad \text{for } i = 1, \dots, n$$
(1)

A text prompt is also mapped into the embedding space, resulting in a sequence of language token embeddings X^l . The LLM takes an input in which visual embeddings are interleaved with language embeddings to generate output Y:

$$Y = \mathcal{M}_L(X_1^v, \dots, X_n^v, X^l) \tag{2}$$



Figure 2: (a) The architecture of M4CXR. Utilizing the LLaVA framework, it allows visual tokens from each image to be inserted at designated positions among the text tokens. (b) Schema for constructing a CXR visual instruction-following dataset. Diverse tasks of three types are combined with appropriate sampling ratios.



Figure 3: Example of multi-turn CoT prompting. M4CXR first identifies findings in the CXR image, then generates a report. Findings and corresponding sentences are colormatched for readability.

The output of the MLLM $Y = \{y_t\}_{t=1}^T$, consisting of T language tokens, is generated auto-regressively:

$$p(Y|X_1^v, \dots, X_n^v, X^l) = \prod_{t=1}^T p(y_t|X_1^v, \dots, X_n^v, X^l, y_{< t}).$$
 (3)

Bounding Box Representation. To enable visual grounding, MLLMs need to represent spatial information as text tokens. Referring to You et al. (2024), we use bounding box coordinates enclosed in square brackets without any additional special tokens. The coordinates $[x_1, y_1, x_2, y_2]$ represent the top-left (x_1, y_1) and bottom-right (x_2, y_2) points of the bounding box on the image. As input images are preprocessed to a uniform size by the vision encoder, the bounding box coordinates are also normalized (Chen et al. 2023a) to integer values between 0 and 100.

Multi-turn Chain-of-Thought Prompting

We derive insights from PromptMRG (Jin et al. 2024) to enhance the clinical accuracy of MRG. It classifies observable

lesions from a CXR image and then uses the result as input prompts to generate reports. We follow a similar two-step process, first identifying key observations and then generating reports. The difference from PromptMRG is that, in our work, a single MLLM performs both classification and report generation sequentially, enabling end-to-end learning.

The gist is to divide the input prompt into two questions, creating multi-turn conversational data. In the first question, we present a list of potential observations and prompt the model to identify those visible in the CXR image. Then, we ask the model to generate a report based on the prior conversation. This approach resembles the reasoning process of a human radiologist, who first identifies lesions visible in the image and then generates a diagnostic report. It can be seen as a variant of the commonly used CoT (Wei et al. 2022) prompting, hence we call it CoT MRG. Figure 3 shows the multi-turn CoT prompting process. To validate the effective-ness of CoT MRG, we conduct a comparative experiment with classification and MRG as separate tasks.

Integrating Chest X-ray Interpretation Tasks

We focus on the following features of an MLLM: the capability for reasoning through multi-turn chat, the ability to utilize multiple images due to its long context length, and the flexibility for multi-task learning. To exploit these features, we construct a CXR visual instruction-following dataset, as schematized in Figure 2 (b). We define a task set across three task types: MRG, image understanding, and VQA, and transform the corresponding datasets into conversational data. Detailed task descriptions and conversation templates are provided in the Appendix (Table 8, 9).

Report Generation: Various Scenarios. A single radiologic study can include multiple CXR images taken from different views: posterior-anterior (PA), anterior-posterior (AP), or lateral. Additionally, when the same patient undergoes follow-up studies, information from prior studies can also be utilized in the radiologic interpretation.

In contrast to previous studies that generate reports from a single image (Jin et al. 2024; Chaves et al. 2024; Hyland et al. 2024; Tanida et al. 2023), we consider MRG tasks in three different scenarios based on the available inputs. The *single-image* scenario incorporates one image from a single study. The *multi-image* scenario accesses multiple images from different views within a single study. The *multi-study* scenario accepts images and the corresponding report from prior study as inputs, along with follow-up study images. These scenarios have different conversation templates. In this work, since we distinguish tasks by templates, they are treated as separate tasks. Through these different scenario tasks, we expect the model to understand and fully utilize the given information, enabling it to generate more accurate reports when additional inputs are available.

Image Understanding. Disease classification is one of the core tasks in CXR image understanding. We use class label data from various classification datasets. Since each dataset has its own set of disease labels, we specify the relevant labels for each question in the prompt. The model is prompted to identify findings in the image from a list of diseases.

Our task set also includes fine-grained tasks requiring detailed analysis of specific CXR regions. These tasks involve grounding and identification of findings, phrases, organs, and anatomical structures, as well as abnormality detection. We convert bounding box labels and clinical text reports into instruction-following data to help the model link CXR images with spatial information for detailed analysis.

Visual Question Answering. To train M4CXR, we convert all the data into an instruction-following format using predefined templates. However, training the LLM solely on templates presents a risk of overfitting: the model might lose its inherent conversational abilities and respond only in a fixed format. Therefore, we also capitalize on VQA datasets which are already in a conversational format, expecting the model to freely respond to a diverse range of questions.

Experiments

Training Datasets

We collect and integrate various datasets according to the tasks described in the **Methods** section.

Report Generation. MIMIC-CXR (Johnson et al. 2019) encompasses a diverse collection of CXR images along with detailed radiology reports. These reports are extensively annotated, facilitating advanced medical image analysis. For MRG tasks, our focus is on generating the FINDINGS section, which offers an in-depth description of significant observations identified in the CXR images.

We employ CheXbert (Smit et al. 2020) to extract observation labels from the FINDINGS section. CheXbert outputs 4-class (positive, negative, uncertain, blank) classification results for 14 predefined observation labels. For binary classification, all non-positive classes are treated as negative. Its label set of observations is provided as the candidate findings in the first question of CoT prompting.

For the three different MRG scenario tasks, image-report pairs are organized as follows. For *single-image*, every image is used as a data instance. In *multi-image*, images from different views that share the same study ID are gathered to compose study-level data. *Multi-study* combines two consecutive studies of a patient, with the studies arranged in chronological order.

Image Understanding. BRAX (Reis et al. 2022) and CheXpert (Irvin et al. 2019) datasets are used for disease classification. We incorporate datasets that contain bounding boxes along with class labels (VinDr-CXR (Nguyen et al. 2022), ChestX-ray14 (Wang et al. 2017), ChestX-Det10 (Liu, Lian, and Yu 2020), JSRT (Shiraishi et al. 2000), SIIM (Zawacki et al. 2019), RSNA (Shih et al. 2019), COVID-19 Radiography (Chowdhury et al. 2020), COVID-QU-Ex (Rahman et al. 2021; Tahir et al. 2021; Chowdhury et al. 2020), QaTa-COV19 (Degerli et al. 2021)) for various finegrained image understanding tasks in addition to disease classification. MS-CXR (Boecking et al. 2022b) is a dataset that includes image-sentence pairs of bounding boxes and corresponding phrases for images and reports from MIMIC-CXR. ImaGenome (Wu et al. 2021), which is also derived from MIMIC-CXR, annotates text in the report with bounding boxes aligned with 29 anatomical regions.

Visual Question Answering. MIMIC-CXR-VQA (Bae et al. 2024) and MIMIC-Diff-VQA (Hu et al. 2023) datasets are designed to address a wide range of questions based on MIMIC-CXR images. These datasets encompass diverse and comprehensive question-answer pairs, enabling models to effectively handle various CXR-related queries. RaDialog (Pellegrini et al. 2023) is a visual instruction-following dataset designed to facilitate tasks that require structured dialogues. It guides models in understanding and responding to instructions based on CXR images.

Test Datasets

The test set used for evaluating all MRG tasks is from MIMIC-CXR. The official test split includes 2,461 images in the frontal (PA, AP) view and 1,397 in other views, totaling 3,858 images. Each study corresponding to the 2,461 frontal images is considered a separate data instance in *multi-image* and *multi-study*, with the prior study included for the latter if available. As a result, the test set size for *multi-image* and *multi-study* is 2,461, while *single-image* considers two test set sizes: 2,461 and 3,858.

For phrase grounding evaluation, the test set of MS-CXR is used. Each query phrase in MS-CXR corresponds to a single bounding box. To ensure a rigorous evaluation of the task, we exclude all MS-CXR images from the training set of other datasets across all tasks. Following the data split of Chen et al. (2023b), the test set has 167 images.

We employ the test set of MIMIC-CXR-VQA for evaluating medical VQA. Of the 13,793 total samples, we exclude 2,484 samples that lack answers, resulting in 11,309 samples for assessment. Additionally, we incorporate SLAKE (Liu et al. 2021), which is not included in the training data. From its test set of 2,094 samples, we use only CXR images with close-ended questions in English, retaining 114 samples.

Evaluation Metrics

Report Generation. The generated reports are assessed using natural language generation (NLG) and clinical met-

Model	CheXbert					Rad-Graph	BL	EU	ROUGE	Test set
in a citati	mF1-14	mF1-5	MF1-14	MF1-5	eF1-14	F1	-1	-4	-L	size
LLM-CXR [†]	36.0	-	21.1	-	-	-	9.2	1.5	16.2	3,530
RaDialog	-	-	39.4	-	-	-	34.6	9.5	27.1	-
METransformer[†]	-	-	-	-	31.1	-	38.6	12.4	29.1	3,269
DCL^\dagger	-	-	-	-	37.3	-	-	10.9	28.4	-
PromptMRG	-	-	38.1	-	47.6	-	39.8	11.2	26.8	3,858
LM-RRG	-	-	-	-	48.4	-	-	12.2	29.6	-
Med-PaLM M 84B	53.6	57.9	39.8	51.6	-	26.7	32.3	11.5	27.5	4,834
CheXagent*	39.3	41.2	24.7	34.5	-	-	16.9	4.7	21.5	2,461
MAIRA-1	55.7	56.0	38.6	47.7	-	24.3	39.2	14.2	28.9	2,461
LLaVA-Rad	57.3	57.4	39.5	47.7	-	-	38.1	15.4	30.6	2,461
M4CXR	58.1 60.6	61.6 61.8	38.8 40.0	49.5 49.5	50.2 53.6	21.7 21.8	33.3 33.9	10.2 10.3	28.4 28.5	3,858 2,461

Table 1: *Single-image* MRG performance on the MIMIC-CXR test set. The best performances are highlighted in bold. The numbers for CheXagent* were obtained from Chaves et al. (2024). Models marked with [†] report classification F1 scores using a different labeler, CheXpert (Irvin et al. 2019).

rics. We use the traditional BLEU (Papineni et al. 2002) and ROUGE-L (Lin 2004) for NLG metrics.

Regarding clinical accuracy, we calculate F1 scores from CheXbert classification results. Macro-averaged F1 (MF1) and micro-averaged F1 (mF1) as well as example-based average F1 (eF1) scores are computed. F1-14 indicates the F1 score calculated over all 14 CheXbert labels, while F1-5 is computed for only 5 of those labels (cardiomegaly, edema, consolidation, atelectasis, pleural effusion).

RadGraph-F1 (Yu et al. 2023) is another clinical accuracy metric. RadGraph (Jain et al. 2021) extracts a graph composed of clinical entities and relations from a radiology report. RadGraph-F1 is the average of F1 scores calculated for both the entities and relations of the graph.

Medical Phrase Grounding. To evaluate the performance of phrase grounding, we calculate the intersection over union (IoU) between the ground-truth and predicted bounding boxes. We use the mean IoU (mIoU) averaged over all data as the evaluation metric. Additionally, accuracy is determined by considering predicted boxes with an IoU of 0.5 or higher as correct predictions.

Visual Question Answering. Accuracy, recall, and BLEU-1 are used to assess the VQA performance. Exact matches between predictions and ground-truth are counted to calculate the accuracy. To consider partial matches for open-ended responses, recall is calculated by measuring the proportion of ground-truth words present in the generated sequences. For MIMIC-CXR-VQA, which includes open-ended questions, BLEU-1 scores are also calculated for each test sample and then averaged.

Model Training

We use a randomly initialized C-Abstractor as the projector to efficiently handle multiple radiology images and reduce image tokens. The LLM and vision encoder are pretrained models: Mistral-7B-Instruct-v0.2 (Jiang et al. 2023) and Rad-DINO (Pérez-García et al. 2024), respectively. The training involves two stages. First, we pre-train the projector while keeping the LLM and vision encoder frozen, using only CXR images and reports without instruction prompts. Subsequently, the vision encoder, projector, and LLM are trained together for visual instruction tuning. We apply LoRA (Hu et al. 2021) for LLM training to reduce computational cost. Each instruction-following data input includes images and instruction texts, and the model is trained to predict the corresponding responses using cross-entropy loss.

The sampling ratio of datasets for each task significantly impacts M4CXR's performance. The empirical ratio for multi-tasking capabilities was determined to be approximately 54%, 35%, and 11% for MRG, image understanding, and VQA, respectively. The detailed training hyperparameters (Table 6), specific sampling ratios (Table 10), and the exploration of sampling ratios are provided in the Appendix.

Results

Medical Report Generation

Single-Image. To evaluate *single-image* MRG performance, we compared our model with state-of-the-art MRG models, including LLM-CXR, RaDialog, METransformer (Wang et al. 2023), DCL (Li et al. 2023b), PromptMRG, LM-RRG (Zhou et al. 2024), CheXagent, MAIRA-1 (Hyland et al. 2024), Med-PaLM M, and LLaVA-Rad (Chaves et al. 2024). The results are summarized in Table 1. For direct comparison with other models, we evaluated M4CXR using two distinct test sets: one with 2,461 frontal images and another with 3,858 images of all views. Including lateral views decreased clinical accuracy, suggesting the complexity of recognizing observations from lateral images.

Among the models evaluated on only frontal images, M4CXR attained the highest CheXbert clinical accuracy, with mF1-14 and MF1-14 scores of 60.6 and 40.0, respectively. In the evaluation that included all views, our model outperformed PromptMRG with an eF1-14 of 50.2. Although Med-PaLM M used a different all-view test set,

which limits direct comparison, our model showed a significant advantage in mF1, with similar or slightly lower MF1s.

However, M4CXR did not achieve the best scores in both RadGraph-F1 and NLG metrics. To improve clinical accuracy, CoT prompting provides candidate findings in the first question before generating reports. This approach may have led the model to use terms from the findings list rather than words that exactly match those in the ground-truth report, which could be a contributing factor to the observed lower NLG scores. RadGraph-F1, calculated by extracting entities and relations from the report, was likely reduced due to this mismatch.

Various Scenarios. Table 2 shows the evaluation results across various input scenarios. *Multi-image* improved clinical accuracy compared to *single-image*, with CheXbert mF1-14 and MF1-14 scores of 61.1 and 41.0, respectively. Providing the prior study as additional input improved MF1, though mF1 slightly decreased. These results demonstrate that our model can effectively utilize available inputs to generate medical reports in various scenarios.

Although test conditions differ, we compared our model's performance with two baseline models. (Tu et al. 2024) reported the zero-shot generalization results of Med-PaLM M in a two-view setting. Since it was trained only on *singleimage*, its effectiveness diminished, unlike M4CXR leveraging additional images for enhanced outcomes. MAIRA-2 accepts a prior study and multiple images as input, akin to our *multi-study* scenario. Its performance matches ours closely: MAIRA-2 achieved a higher MF1-14 of 42.7, whereas M4CXR surpassed it with an mF1-14 of 60.7.

Model	MRG Scenario	CheXbert		
widder	WIKO Sechario	mF1-14	MF1-14	
Med-PaLM M 84B	Multi-image	50.5	37.8	
MAIRA-27B	Multi-study	58.5	42.7	
	Single-image	60.6	40.0	
M4CXR	Multi-image	61.1	41.0	
	Multi-study	60.7	42.0	

Table 2: Evaluation results of MRG performance in various scenarios.

Medical Phrase Grounding

Table 3 compares the medical phrase grounding results of TransVG (Deng et al. 2022), MedRPG, MAIRA-2, and M4CXR, on the MS-CXR test set. While MedRPG achieved the highest accuracy and mIoU, MedRPG and TransVG are specialized for phrase grounding and cannot perform other tasks. Given this limitation, M4CXR, a multi-tasking model, shows competitive performance with an accuracy of 68.3 and mIoU of 57.9. MAIRA-2, utilizing a private dataset for grounded report generation, reported a similar mIoU to ours.

Visual Question Answering

Table 4 presents the results of medical VQA evaluation. We conducted our own evaluations on open-source medical MLLMs: LLaVA version of RaDialog (Pellegrini et al.

Model	Accuracy	mIOU
TransVG	65.9	58.9
MedRPG	69.9	59.4
MAIRA-27B	-	57.8
M4CXR	68.3	57.9

Table 3: Referring expression grounding results. Note that MAIRA-2 7B was evaluated on a different test split.

2024), RadFM (Wu et al. 2023), and CheXagent. M4CXR outperformed the other models, except for the recall on MIMIC-CXR-VQA. CheXagent included the dataset in its training, likely explaining its high recall of 72.8. The lower accuracy on MIMIC-CXR-VQA, compared to SLAKE, is due to its greater variety of questions. Moreover, the higher BLEU-1 score compared to accuracy suggests the presence of partially correct answers.

Model	MI	MIC-CX	SLAKE		
11100001	Acc	Recall	BLEU-1	Acc	Recall
RaDialog	0.0	43.0	5.0	0.0	45.6
RadFM	11.2	36.7	13.8	68.4	69.7
CheXagent	59.0	72.8	62.5	71.1	73.2
M4CXR	62.3	70.3	66.4	85.1	86.0

Table 4: Medical VQA performance.

Ablation Study

We examined the effects of multi-turn CoT prompting and the combination of task types through comparative experiments (Table 5). Generating reports without CoT prompting in M4CXR led to a significant drop in clinical accuracy. Exp1, trained on single-turn MRG without the first question of CoT prompts, yielded similar results. Exp2 involved training on single-turn MRG and disease classification as separate tasks. Testing Exp2 with CoT prompts improved clinical accuracy, but NLG metrics fell significantly, with the generated text being merely a list of identified observations. The experiments reveal the effectiveness of multi-turn CoT prompting in MRG's training and testing phases.

Exp3 and Exp4 excluded image understanding and VQA tasks, respectively. In Exp3, while MRG performance was on par and VQA performance improved due to the increased proportion of VQA data, the absence of image understanding data precluded visual grounding. In Exp4, the BLEU-1 score for VQA dropped significantly from 66.4 to 36.1, suggesting the necessity of VQA datasets in understanding and answering free-form questions. Moreover, the slight decrease in Exp4's clinical accuracy compared to M4CXR implies the contribution of VQA datasets to MRG performance improvement. Detailed examples are in the Appendix.

Qualitative Analysis

Figure 4 shows two examples of dialogs between a user and M4CXR, each for visual grounding and VQA. Additional

		Training		Test	Che	Xbert	N	ILG	Grounding	VQA
	MRG	ImgUnd	VQA	CoT	mF1-14	MF1-14	BLEU-4	ROUGE-L	mIoU	BLEU-1
M4CXR	1	1	1	✓ ×	60.6 49.2	40.0 31.2	10.3 9.3	28.5 27.3	57.9	66.4
Exp1	\checkmark^{s}	1	1	×	50.0	33.1	9.3	27.7	56.1	66.1
Exp2	$\checkmark^{s,c}$	1	1	✓ ×	58.5 49.2	37.5 31.6	0.6 9.2	10.9 27.2	56.3	66.6
Exp3	1		1	1	60.6	38.2	10.1	28.2	-	67.9
Exp4	1	1		1	59.3	38.3	10.1	28.0	57.9	36.1

Table 5: Results of ablation studies. "ImgUnd" stands for "Image Understanding". \checkmark^s and $\checkmark^{s,c}$ indicate single-turn MRG without CoT prompting and the training of single-turn MRG and disease classification as separate tasks, respectively.



Figure 4: Examples of M4CXR's performance in (a) visual grounding and (b) VQA. The images are selected from the test splits of MS-CXR and MIMIC-CXR, respectively.

examples can be found in the Appendix. In Figure 4 (a), after the multi-turn CoT prompting to generate a report, the user asks a third question to locate the region of a sentence in the generated report. The model provides the coordinates of a bounding box, depicted as a red dashed box. The prediction covers all of the ground-truth, represented by the yellow solid box, possessing sufficient explanatory power to indicate what the corresponding sentence represents.

Figure 4 (b) illustrates M4CXR's ability to write a clear and concise summary from the generated report. On the left, the ground-truth report is provided for comparison. For enhanced visual clarity, findings and corresponding sentences are color-matched. M4CXR accurately identifies all the observations in the first question, and subsequently generates a comprehensive report without omitting any detail. In response to the third question which asks for a one-sentence summary, the model offers a concise statement that includes every noted finding. A notable flaw is that expressions such as "increasing" and "worsening," which would be meaningful only in a *multi-study* scenario, were generated despite the *single-image* context. While it is understandable that this hallucination is due to the ground-truth reports containing such expressions during training, it is evident that there is room for improvement.

Conclusion

In this study, we introduced M4CXR, a multi-modal LLM aimed at enhancing CXR interpretation by leveraging the versatile advantages of LLMs. Trained on a CXR visual instruction-following dataset constructed by appropriately combining various datasets, M4CXR is capable of performing multiple tasks. Exploiting its reasoning capability, we proposed a novel chain-of-thought prompting strategy that significantly improved the clinical accuracy of report generation. By identifying observations from X-ray images and generating descriptions through multi-turn conversations, M4CXR achieved notable improvements in CheXbert F1 scores compared to existing models. The model's ability to handle multiple images and incorporate prior studies demonstrated its adaptability in diverse clinical scenarios.

Beyond report generation, M4CXR also proved to be highly effective in visual grounding and VQA. M4CXR competed well with models specialized for visual grounding, while also exhibiting outstanding performance in VQA. Qualitative analysis highlights M4CXR's flexibility in answering free-form questions. With future work to address limitations such as low NLG metrics and the presence of hallucinations, we anticipate that further advancements in M4CXR will contribute to the development of a highly capable radiology assistant.

References

Alayrac, J.-B.; Donahue, J.; Luc, P.; Miech, A.; Barr, I.; Hasson, Y.; Lenc, K.; Mensch, A.; Millican, K.; Reynolds, M.; et al. 2022. Flamingo: a visual language model for few-shot learning. *Advances in neural information processing systems*, 35: 23716–23736.

Bae, S.; Kyung, D.; Ryu, J.; Cho, E.; Lee, G.; Kweon, S.; Oh, J.; Ji, L.; Chang, E.; Kim, T.; et al. 2024. EHRXQA: A multi-modal question answering dataset for electronic health records with chest x-ray images. *Advances in Neural Information Processing Systems*, 36.

Bannur, S.; Bouzid, K.; Castro, D. C.; Schwaighofer, A.; Bond-Taylor, S.; Ilse, M.; Pérez-García, F.; Salvatelli, V.; Sharma, H.; Meissen, F.; Ranjit, M.; Srivastav, S.; Gong, J.; Falck, F.; Oktay, O.; Thieme, A.; Lungren, M. P.; Wetscherek, M. T.; Alvarez-Valle, J.; and Hyland, S. L. 2024. MAIRA-2: Grounded Radiology Report Generation. arXiv:2406.04449.

Boecking, B.; Usuyama, N.; Bannur, S.; Castro, D. C.; Schwaighofer, A.; Hyland, S.; Wetscherek, M.; Naumann, T.; Nori, A.; Alvarez-Valle, J.; et al. 2022a. Making the most of text semantics to improve biomedical vision–language processing. In *European conference on computer vision*, 1– 21. Springer.

Boecking, B.; Usuyama, N.; Bannur, S.; Coelho de Castro, D.; Schwaighofer, A.; Hyland, S.; Wetscherek, M. T.; Naumann, T.; Nori, A.; Alvarez Valle, J.; Poon, H.; and Oktay, O. 2022b. MS-CXR: Making the Most of Text Semantics to Improve Biomedical Vision-Language Processing (version 0.1). https://doi.org/10.13026/b90j-vb87.

Cha, J.; Kang, W.; Mun, J.; and Roh, B. 2024. Honeybee: Locality-enhanced projector for multimodal llm. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 13817–13827.

Chaves, J. M. Z.; Huang, S.-C.; Xu, Y.; Xu, H.; Usuyama, N.; Zhang, S.; Wang, F.; Xie, Y.; Khademi, M.; Yang, Z.; Awadalla, H.; Gong, J.; Hu, H.; Yang, J.; Li, C.; Gao, J.; Gu, Y.; Wong, C.; Wei, M.; Naumann, T.; Chen, M.; Lungren, M. P.; Chaudhari, A.; Yeung-Levy, S.; Langlotz, C. P.; Wang, S.; and Poon, H. 2024. Towards a clinically accessible radiology foundation model: open-access and lightweight, with automated evaluation. arXiv:2403.08002.

Chen, J.; Zhu, D.; Shen, X.; Li, X.; Liu, Z.; Zhang, P.; Krishnamoorthi, R.; Chandra, V.; Xiong, Y.; and Elhoseiny, M. 2023a. MiniGPT-v2: large language model as a unified interface for vision-language multi-task learning. arXiv:2310.09478.

Chen, Z.; Song, Y.; Chang, T.-H.; and Wan, X. 2022. Generating Radiology Reports via Memory-driven Transformer. arXiv:2010.16056.

Chen, Z.; Varma, M.; Delbrouck, J.-B.; Paschali, M.; Blankemeier, L.; Veen, D. V.; Valanarasu, J. M. J.; Youssef, A.; Cohen, J. P.; Reis, E. P.; Tsai, E. B.; Johnston, A.; Olsen, C.; Abraham, T. M.; Gatidis, S.; Chaudhari, A. S.; and Langlotz, C. 2024. CheXagent: Towards a Foundation Model for Chest X-Ray Interpretation. arXiv:2401.12208.

Chen, Z.; Zhou, Y.; Tran, A.; Zhao, J.; Wan, L.; Ooi, G. S. K.; Cheng, L. T.-E.; Thng, C. H.; Xu, X.; Liu, Y.; et al. 2023b. Medical phrase grounding with region-phrase context contrastive alignment. In *International Conference on Medical Image Computing and Computer-Assisted Intervention*, 371–381. Springer.

Chowdhury, M. E. H.; Rahman, T.; Khandakar, A.; Mazhar, R.; Kadir, M. A.; Mahbub, Z. B.; Islam, K. R.; Khan, M. S.; Iqbal, A.; Emadi, N. A.; Reaz, M. B. I.; and Islam, M. T. 2020. Can AI Help in Screening Viral and COVID-19 Pneumonia? *IEEE Access*, 8: 132665–132676.

Dao, T. 2023. FlashAttention-2: Faster Attention with Better Parallelism and Work Partitioning. arXiv:2307.08691.

Degerli, A.; Ahishali, M.; Yamac, M.; Kiranyaz, S.; Chowdhury, M. E. H.; Hameed, K.; Hamid, T.; Mazhar, R.; and Gabbouj, M. 2021. COVID-19 infection map generation and detection from chest X-ray images. *Health Information Science and Systems*, 9(1): 15.

Deng, J.; Yang, Z.; Chen, T.; Zhou, W.; and Li, H. 2022. TransVG: End-to-End Visual Grounding with Transformers. arXiv:2104.08541.

Hu, E. J.; Shen, Y.; Wallis, P.; Allen-Zhu, Z.; Li, Y.; Wang, S.; Wang, L.; and Chen, W. 2021. LoRA: Low-Rank Adaptation of Large Language Models. arXiv:2106.09685.

Hu, X.; Gu, L.; An, Q.; Zhang, M.; Liu, L.; Kobayashi, K.; Harada, T.; Summers, R. M.; and Zhu, Y. 2023. Expert Knowledge-Aware Image Difference Graph Representation Learning for Difference-Aware Medical Visual Question Answering. In *Proceedings of the 29th ACM SIGKDD Conference on Knowledge Discovery and Data Mining*, KDD '23. ACM.

Hyland, S. L.; Bannur, S.; Bouzid, K.; Castro, D. C.; Ranjit, M.; Schwaighofer, A.; Pérez-García, F.; Salvatelli, V.; Srivastav, S.; Thieme, A.; Codella, N.; Lungren, M. P.; Wetscherek, M. T.; Oktay, O.; and Alvarez-Valle, J. 2024. MAIRA-1: A specialised large multimodal model for radiology report generation. arXiv:2311.13668.

Irvin, J.; Rajpurkar, P.; Ko, M.; Yu, Y.; Ciurea-Ilcus, S.; Chute, C.; Marklund, H.; Haghgoo, B.; Ball, R.; Shpanskaya, K.; Seekins, J.; Mong, D. A.; Halabi, S. S.; Sandberg, J. K.; Jones, R.; Larson, D. B.; Langlotz, C. P.; Patel, B. N.; Lungren, M. P.; and Ng, A. Y. 2019. CheXpert: A Large Chest Radiograph Dataset with Uncertainty Labels and Expert Comparison. arXiv:1901.07031.

Jain, S.; Agrawal, A.; Saporta, A.; Truong, S.; Duong, D. N. D. N.; Bui, T.; Chambon, P.; Zhang, Y.; Lungren, M.; Ng, A.; Langlotz, C.; Rajpurkar, P.; and Rajpurkar, P. 2021. RadGraph: Extracting Clinical Entities and Relations from Radiology Reports. In Vanschoren, J.; and Yeung, S., eds., *Proceedings of the Neural Information Processing Systems Track on Datasets and Benchmarks*, volume 1. Jiang, A. Q.; Sablayrolles, A.; Mensch, A.; Bamford, C.; Chaplot, D. S.; de las Casas, D.; Bressand, F.; Lengyel, G.; Lample, G.; Saulnier, L.; Lavaud, L. R.; Lachaux, M.-A.; Stock, P.; Scao, T. L.; Lavril, T.; Wang, T.; Lacroix, T.; and Sayed, W. E. 2023. Mistral 7B. arXiv:2310.06825.

Jin, H.; Che, H.; Lin, Y.; and Chen, H. 2024. PromptMRG: Diagnosis-Driven Prompts for Medical Report Generation. *Proceedings of the AAAI Conference on Artificial Intelligence*, 38(3): 2607–2615.

Jing, B.; Xie, P.; and Xing, E. 2018. On the Automatic Generation of Medical Imaging Reports. In *Proceedings of the* 56th Annual Meeting of the Association for Computational Linguistics (Volume 1: Long Papers). Association for Computational Linguistics.

Johnson, A. E. W.; Pollard, T. J.; Berkowitz, S. J.; Greenbaum, N. R.; Lungren, M. P.; Deng, C.; Mark, R. G.; and Horng, S. 2019. MIMIC-CXR: A large publicly available database of labeled chest radiographs. *CoRR*, abs/1901.07042.

Lee, S.; Kim, W. J.; Chang, J.; and Ye, J. C. 2024. LLM-CXR: Instruction-Finetuned LLM for CXR Image Understanding and Generation. arXiv:2305.11490.

Li, J.; Li, D.; Savarese, S.; and Hoi, S. 2023a. Blip-2: Bootstrapping language-image pre-training with frozen image encoders and large language models. In *International conference on machine learning*, 19730–19742. PMLR.

Li, M.; Lin, B.; Chen, Z.; Lin, H.; Liang, X.; and Chang, X. 2023b. Dynamic graph enhanced contrastive learning for chest x-ray report generation. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 3334–3343.

Lin, C.-Y. 2004. Rouge: A package for automatic evaluation of summaries. In *Text summarization branches out*, 74–81.

Liu, B.; Zhan, L.; Xu, L.; Ma, L.; Yang, Y.; and Wu, X. 2021. SLAKE: A Semantically-Labeled Knowledge-Enhanced Dataset for Medical Visual Question Answering. *CoRR*, abs/2102.09542.

Liu, H.; Li, C.; Wu, Q.; and Lee, Y. J. 2024. Visual instruction tuning. *Advances in neural information processing systems*, 36.

Liu, J.; Lian, J.; and Yu, Y. 2020. ChestX-Det10: Chest X-ray Dataset on Detection of Thoracic Abnormalities. arXiv:2006.10550.

Nguyen, H. Q.; Lam, K.; Le, L. T.; Pham, H. H.; Tran, D. Q.; Nguyen, D. B.; Le, D. D.; Pham, C. M.; Tong, H. T. T.; Dinh, D. H.; Do, C. D.; Doan, L. T.; Nguyen, C. N.; Nguyen, B. T.; Nguyen, Q. V.; Hoang, A. D.; Phan, H. N.; Nguyen, A. T.; Ho, P. H.; Ngo, D. T.; Nguyen, N. T.; Nguyen, N. T.; Dao, M.; and Vu, V. 2022. VinDr-CXR: An open dataset of chest X-rays with radiologist's annotations. arXiv:2012.15029.

Papineni, K.; Roukos, S.; Ward, T.; and Zhu, W.-J. 2002. Bleu: a method for automatic evaluation of machine translation. In *Proceedings of the 40th annual meeting of the Association for Computational Linguistics*, 311–318.

Pellegrini, C.; Özsoy, E.; Busam, B.; Navab, N.; and Keicher, M. 2023. RaDialog: A Large Vision-Language Model for Radiology Report Generation and Conversational Assistance. arXiv:2311.18681.

Pellegrini, C.; Özsoy, E.; Busam, B.; Navab, N.; and Keicher, M. 2024. RaDialog Instruct Dataset (version 1.1.0). PhysioNet.

Pérez-García, F.; Sharma, H.; Bond-Taylor, S.; Bouzid, K.; Salvatelli, V.; Ilse, M.; Bannur, S.; Castro, D. C.; Schwaighofer, A.; Lungren, M. P.; Wetscherek, M.; Codella, N.; Hyland, S. L.; Alvarez-Valle, J.; and Oktay, O. 2024. RAD-DINO: Exploring Scalable Medical Image Encoders Beyond Text Supervision. arXiv:2401.10815.

Raffel, C.; Shazeer, N.; Roberts, A.; Lee, K.; Narang, S.; Matena, M.; Zhou, Y.; Li, W.; and Liu, P. J. 2020. Exploring the limits of transfer learning with a unified text-to-text transformer. *Journal of machine learning research*, 21(140): 1–67.

Rahman, T.; Khandakar, A.; Qiblawey, Y.; Tahir, A.; Kiranyaz, S.; Abul Kashem, S.; Islam, M.; Al Maadeed, S.; Zughaier, S.; Khan, M.; and Chowdhury, M. 2021. Exploring the Effect of Image Enhancement Techniques on COVID-19 Detection using Chest X-rays Images. *Computers in Biology and Medicine*, 104319.

Reis, E. P.; de Paiva, J. P.; da Silva, M. C.; Ribeiro, G. A.; Paiva, V. F.; Bulgarelli, L.; Lee, H. M.; Santos, P. V.; Brito, V. M.; Amaral, L. T.; Beraldo, G. L.; Haidar Filho, J. N.; Teles, G. B.; Szarf, G.; Pollard, T.; Johnson, A. E.; Celi, L. A.; and Amaro, E. J. 2022. BRAX, Brazilian labeled chest x-ray dataset. *Scientific Data*, 9(1): 487.

Shih, G.; Wu, C.; Halabi, S.; Kohli, M.; Prevedello, L.; Cook, T.; Sharma, A.; Amorosa, J.; Arteaga, V.; Galperin-Aizenberg, M.; Gill, R.; Godoy, M.; Hobbs, S.; Jeudy, J.; Laroia, A.; Shah, P.; Vummidi, D.; Yaddanapudi, K.; and Stein, A. 2019. Augmenting the national institutes of health chest radiograph dataset with expert annotations of possible pneumonia. *Radiology: Artificial Intelligence*, 1(1).

Shin, H.-C.; Roberts, K.; Lu, L.; Demner-Fushman, D.; Yao, J.; and Summers, R. M. 2016. Learning to Read Chest X-Rays: Recurrent Neural Cascade Model for Automated Image Annotation. In *Proceedings of the IEEE Conference on Computer Vision and Pattern Recognition (CVPR)*.

Shiraishi, J.; Katsuragawa, S.; Ikezoe, J.; Matsumoto, T.; Kobayashi, T.; Komatsu, K.; Matsui, M.; Fujita, H.; Kodera, Y.; and Doi, K. 2000. Development of a digital image database for chest radiographs with and without a lung nodule: receiver operating characteristic analysis of radiologists' detection of pulmonary nodules. *AJR Am J Roentgenol*, 174(1): 71–74.

Smit, A.; Jain, S.; Rajpurkar, P.; Pareek, A.; Ng, A. Y.; and Lungren, M. P. 2020. CheXbert: Combining Automatic Labelers and Expert Annotations for Accurate Radiology Report Labeling Using BERT. arXiv:2004.09167.

Tahir, A.; Chowdhury, M.; Qiblawey, Y.; Khandakar, A.; Rahman, T.; Kiranyaz, S.; Khurshid, U.; Ibtehaz, N.; Mahmud, S.; and Ezeddin, M. 2021. COVID-QU-Ex Dataset. https://www.kaggle.com/datasets/tahirahmed/covidquex.

Tanida, T.; Müller, P.; Kaissis, G.; and Rueckert, D. 2023. Interactive and explainable region-guided radiology report generation. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 7433–7442.

Thawkar, O.; Shaker, A.; Mullappilly, S. S.; Cholakkal, H.; Anwer, R. M.; Khan, S.; Laaksonen, J.; and Khan, F. S. 2023. XrayGPT: Chest Radiographs Summarization using Medical Vision-Language Models. arXiv:2306.07971.

Tu, T.; Azizi, S.; Driess, D.; Schaekermann, M.; Amin, M.; Chang, P.-C.; Carroll, A.; Lau, C.; Tanno, R.; Ktena, I.; et al. 2024. Towards generalist biomedical AI. *NEJM AI*, 1(3): AIoa2300138.

Wang, X.; Peng, Y.; Lu, L.; Lu, Z.; Bagheri, M.; and Summers, R. M. 2017. ChestX-Ray8: Hospital-Scale Chest X-Ray Database and Benchmarks on Weakly-Supervised Classification and Localization of Common Thorax Diseases. In 2017 IEEE Conference on Computer Vision and Pattern Recognition (CVPR), 3462–3471.

Wang, Z.; Liu, L.; Wang, L.; and Zhou, L. 2023. Metransformer: Radiology report generation by transformer with multiple learnable expert tokens. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 11558–11567.

Wei, J.; Wang, X.; Schuurmans, D.; Bosma, M.; Xia, F.; Chi, E.; Le, Q. V.; Zhou, D.; et al. 2022. Chain-of-thought prompting elicits reasoning in large language models. *Advances in neural information processing systems*, 35: 24824–24837.

Wu, C.; Zhang, X.; Zhang, Y.; Wang, Y.; and Xie, W. 2023. Towards Generalist Foundation Model for Radiology by Leveraging Web-scale 2D&3D Medical Data. arXiv:2308.02463.

Wu, J. T.; Agu, N. N.; Lourentzou, I.; Sharma, A.; Paguio, J. A.; Yao, J. S.; Dee, E. C.; Mitchell, W.; Kashyap, S.; Giovannini, A.; Celi, L. A.; and Moradi, M. 2021. Chest ImaGenome Dataset for Clinical Reasoning. arXiv:2108.00316.

Yang, L.; Xu, S.; Sellergren, A.; Kohlberger, T.; Zhou, Y.; Ktena, I.; Kiraly, A.; Ahmed, F.; Hormozdiari, F.; Jaroensri, T.; et al. 2024. Advancing multimodal medical capabilities of Gemini. arXiv:2405.03162.

You, H.; Zhang, H.; Gan, Z.; Du, X.; Zhang, B.; Wang, Z.; Cao, L.; Chang, S.-F.; and Yang, Y. 2024. Ferret: Refer and Ground Anything Anywhere at Any Granularity. In *The Twelfth International Conference on Learning Representations*.

You, K.; Gu, J.; Ham, J.; Park, B.; Kim, J.; Hong, E. K.; Baek, W.; and Roh, B. 2023. CXR-CLIP: Toward Large Scale Chest X-ray Language-Image Pre-training. In *Medical Image Computing and Computer Assisted Intervention – MICCAI 2023*, 101–111. Springer Nature Switzerland.

Yu, F.; Endo, M.; Krishnan, R.; Pan, I.; Tsai, A.; Reis, E. P.; Fonseca, E. K. U. N.; Lee, H. M. H.; Abad, Z. S. H.; Ng, A. Y.; et al. 2023. Evaluating progress in automatic chest x-ray radiology report generation. *Patterns*, 4(9).

Zawacki, A.; Wu, C.; Shih, G.; Elliott, J.; Fomitchev, M.; Hussain, M.; ParasLakhani; Culliton, P.; and Bao, S. 2019. SIIM-ACR Pneumothorax Segmentation. Zhang, Y.; Wang, X.; Xu, Z.; Yu, Q.; Yuille, A.; and Xu, D. 2020. When Radiology Report Generation Meets Knowledge Graph. *Proceedings of the AAAI Conference on Artificial Intelligence*, 34(07): 12910–12917.

Zheng, Y.; Gan, W.; Chen, Z.; Qi, Z.; Liang, Q.; and Yu, P. S. 2024. Large Language Models for Medicine: A Survey. arXiv:2405.13055.

Zhou, Z.; Shi, M.; Wei, M.; Alabi, O.; Yue, Z.; and Vercauteren, T. 2024. Large Model driven Radiology Report Generation with Clinical Quality Reinforcement Learning. arXiv:2403.06728.

Zhu, D.; Chen, J.; Shen, X.; Li, X.; and Elhoseiny, M. 2023. MiniGPT-4: Enhancing Vision-Language Understanding with Advanced Large Language Models. arXiv:2304.10592.

Appendix

Implementation Details

The detailed hyperparameters for training are summarized in Table 6. C-Abstractor is configured with a depth of 3, an MLP depth of 2, and a hidden size of 1024, to be employed as the projector. The number of visual tokens is set to 361, corresponding to a 19×19 grid. Accordingly, the output of the vision encoder, which is a sequence of 1369 (37×37) tokens processed by Rad-DINO, is compressed to a length of 361. This reduction in the number of tokens, by approximately one-fourth, enables the LLM to handle a greater number of images more efficiently. In the training of the LLM, we configure the LoRA hyperparameters as follows: rank is set to 8, alpha to 32, and dropout to 0.05. Experiments are conducted on two H100 GPUs, using bfloat16 with automatic mixed precision and flash-attention v2 (Dao 2023). For evaluation, sentences are generated with greedy search decoding. One exception is the qualitative evaluation of VQA: to generate more natural responses, multinomial sampling is applied.

Configuration	Pre-training	Instruction tuning		
Training module	Projector	Vision Encoder, Projector, LLM		
Training steps	2k	40k		
Max steps	10k	100k		
Warmup steps	500			
Global batch size	256	64		
Optimizer	AdamW			
Optimizer hyperparameter	$\beta_1 = 0.9, \beta_2$	$= 0.999, \epsilon = 1e^{-8}$		
Learning rate	3e-4	1e-4		
Learning rate schedule	Cos	ine decay		
Weight decay	1e-2	1e-4		
Gradient clipping	1.0			

Table 6: Training hyperparameters.

Details on Chest X-ray Interpretation Tasks

Detailed descriptions for all the CXR interpretation tasks considered in this work are listed in Table 8. Table 9 provides the conversation templates for each task. Although not shown in the table, a system prompt, which specifies that the AI medical assistant should give helpful and detailed answers, is inserted before the first question in each template. Table 10 lists the datasets used for each task, along with the number of training instances and the sampling ratios applied during training.

Dataset Preprocessing

We follow the official split for every dataset. The additional preprocessing steps, taken for each dataset as required, are described in detail. **MIMIC-CXR.** We extract the FINDINGS section using the official preprocessing code¹. Then, we use the text preprocessing code from Chen et al. (2022) with slight modification to remove special characters, numbering, and extra spaces. We exclude reports that lack a FINDINGS section or have a length of less than 5 characters. For *multi-image*, we consider only studies with at most 5 images. For *multi-study*, we include only cases where the combined number of images in the prior and follow-up studies is at most 10. Out of 2,461 follow-up studies in the test set, only 2,333 studies have a prior study and the remaining 128 studies have no preceding study.

VinDR-CXR. Due to the presence of many similar bounding box labels, overlapping boxes in an image are merged into the smallest possible bounding box if they overlap by more than 50 percent.

JSRT. Circular labels are converted to the smallest possible bounding box labels.

Mask Labeled Datasets. The mask labels are converted to bounding boxes that minimally enclose each mask. This conversion is applied for ChestX-Det10, SIIM, COVID-19 Radiography, COVID-QU-EX, and QaTa-COV19.

COVID-19 Radiography. Images with more than three mask regions are excluded.

QaTa-COV19. Data that overlaps with COVID-QU-EX are removed.

RSNA. Only 'lung opacity' and 'normal' labels are used. According to You et al. (2023), 'lung opacity' is referred to as 'pneumonia'.

ImaGenome. We utilize the preprocessing code from Tanida et al. (2023). As a result, we obtain report text, bounding boxes, and anatomical names for 29 anatomical regions.

RaDialog. We use only the images from MIMIC-CXR. Among the tasks introduced in Pellegrini et al. (2023), we exclude RG (report generation). We utilize version 1.1.0 (Pellegrini et al. 2024) of the dataset.

Exploration of Sampling Ratios

In multi-task learning, the data sampling strategy plays a crucial role, with numerous such strategies being applicable. We experiment with various sampling strategies, as shown in Table 7. The tasks are grouped into three task types, and a dataset used within a single task is referred to as a task-dataset. The per-task-dataset (D1) strategy samples uniformly across all task-datasets. The per-size (D2) strategy samples proportionally based on the size of each task-dataset. After setting equal proportions for each task type (per-task-type), within each task-type, D3 samples per-task-dataset while D4 samples per-size.

In D1, which follows the per-task-dataset strategy, MRG accounts for only 3 out of 50 total task-datasets. This very

¹https://github.com/MIT-LCP/mimic-cxr/tree/master/txt

low sampling ratio for MRG explains its low MRG performance. In contrast, in the per-task-type settings of D3 and D4, one-third of the overall tasks are allocated to MRG, leading to improved MRG performance compared to D1. For grounding performance, D2 and D4 showed inferior results, likely because the training split of MS-CXR, which is the dataset used for evaluation, was barely sampled when using the per-size setting. As for VQA, since all datasets are large in size, per-size sampling (D2) exhibited a high BLEU-1 score. Similar to MRG, VQA also constitutes a small proportion of the total task-dataset (4 out of 50), resulting in lower performance for D1, while D3 and D4 achieved higher scores as per-task-type allocates one-third to VQA.

We aimed to find the appropriate sampling ratios where M4CXR remains clinically accurate, supports visual grounding, and retains the conversational abilities of the LLM. To improve MRG performance, we reduced the ratio allocated to VQA and increased the proportion dedicated to MRG. For image understanding, we began with the persize ratio and made further adjustments to enhance grounding performance. Additionally, to improve the instructionfollowing capabilities of the model, we increased the ratio of RaDialog compared to other VQA datasets. While this adjustment led to a decrease in quantitative VQA performance metrics, it contributed in preserving the LLM's ability to understand and respond to instructions in conversation.

Qualitative Examples

Comparison of MRG Scenarios. Figure 5 presents the results of report generation for the same study across various scenarios, highlighting how clinical accuracy improves as more information becomes available. The ground-truth targets, including atelectasis, pleural effusion, and support devices, are each marked in different colors. In *single-image*, one lateral image is used; in *multi-image*, both AP and lateral images are used; and in *multi-study*, three prior images and the corresponding report are included.

In *single-image*, only the support devices are correctly identified, suggesting that a single lateral image may not provide sufficient information for accurate CXR interpretation. In *multi-image*, the addition of a frontal image enables the correct identification of atelectasis. In *multi-study*, all three targets are successfully identified, possibly due to the prior report already containing these three findings. This demonstrates that M4CXR is capable of understanding and utilizing the information provided in various MRG scenarios.

Visual Grounding. Figure 6 presents additional examples of visual grounding. In both (a) and (b), the user asks a third question after MRG to identify the location referenced by a specific phrase, and M4CXR responds with the coordinates of a bounding box.

In Figure 6 (a), the model predicts the presence of lung opacity and atelectasis. The report contains the phrase "volume loss consistent with right upper lobe collapse", which suggests atelectasis. When asked to ground the sentence containing this phrase, the model identifies an area that largely overlaps with the ground-truth bounding box. In Fig-

ure 6 (b), the generated report explicitly mentions the presence of pneumothorax, and the model correctly grounds the upper part of the right lung. These examples demonstrate that M4CXR can be utilized both for generating reports and for determining the locations referenced by the generated text.

Impact of VQA Datasets. Figures 7, 8, and 9 collectively illustrate the differences in VQA performance between M4CXR and the Exp4 model from **Ablation Study**, emphasizing the impact of incorporating VQA datasets. The generated responses were compared using images selected from the MIMIC-CXR test set.

In Figure 7, when requested to rephrase the generated reports in easy language, M4CXR effectively translates medical jargon into general language. In contrast, Exp4 model merely summarizes the report content, continuing to use medical terminology without making it easier to understand.

Figure 8 presents the comparison results for the task of report summarization. M4CXR successfully condenses the report into a concise, single sentence as requested, notably including all important observations such as cardiomegaly, pleural effusion, and support devices. In contrast, the model trained without VQA datasets oversimplifies the content, mentioning only pleural effusion and omitting other critical details.

Figure 9 shows the responses when recommendations for diagnosis and treatment are requested. M4CXR suggests continued observation and follow-up imaging studies based on the findings, and also recommends treatments such as pleural drainage or chest tube placement if the condition worsens. In contrast, Exp4 model provides a simpler response, suggesting evaluation with CT and the need for monitoring.

These comparisons between the two models demonstrate that VQA datasets play a significant role in maintaining the LLM's ability to follow instructions and provide appropriate responses.

Hallucination. In Figure 9, the medical report generated by M4CXR contains comparative contents such as "In comparison with study" and "decreased", even though it was generated in a *single-image* scenario. Similar issues can be observed in other examples, where the model often refers to comparisons or mentions images that are not provided. During the training of *single-image* MRG, the use of certain ground-truth reports, those written specifically in a *multi-study* context, is likely to have induced these hallucinations. To address this issue, ground-truth reports tailored to each MRG scenario are needed, which could be explored in future work.

	Sampling Strategy	CheXbert		NLG		Grounding	VQA
	Sampling Strategy	mF1-14	MF1-14	BLEU-4	ROUGE-L	mIoU	BLEU-1
M4CXR	empirical ratio	60.6	40.0	10.3	28.5	57.9	66.4
D1	per-task-dataset	52.7	32.1	9.1	26.8	58.5	72.8
D2	per-size	58.5	36.1	9.0	27.2	47.9	77.4
D3	per-task-type & per-task-dataset	59.0	37.1	10.0	28.2	57.9	78.3
D4	per-task-type & per-size	59.1	37.7	9.1	27.3	46.4	78.0

Table 7: Performance comparison between different data sampling strategies.

Task Type	Task	Description				
	Single-image	Generates the FINDINGS section of a report from a single CXR image.				
Medical Report	Multi-image	Generates the FINDINGS section of a report from one or more images, in- cluding images from different views within the same study.				
Generation	Multi-study	Generates the FINDINGS section of a report from one or more studies, com bining current and previous studies from a patient.				
	Disease Classification	Identifies diseases within a single image based on the labels covered by each dataset.				
	Finding Grounding	Provides the bounding box coordinates for a given finding, if detected.				
	Grounded Finding	Identifies the finding corresponding to a given bounding box.				
	Abnormality Detection	Identifies and localizes abnormal regions by providing the corresponding bounding box coordinates.				
	Multi Finding Grounding	Identifies the presence of candidate findings and provides the bounding boucoordinates for each detected finding.				
Image	Organ Grounding	Provides the bounding box coordinates for a given organ.				
Understanding	Grounded Organ	Identifies the organ corresponding to a given bounding box.				
	Grounded Phrase Generation	Generates a radiology report phrase for the region corresponding to a given bounding box.				
	Phrase Grounding	Provides the bounding box coordinates for the region described by a given phrase.				
	Anatomical Region Grounding	Provides the bounding box coordinates for a given anatomical region.				
	Grounded Anatomical Region	Identifies the anatomical region corresponding to a given bounding box.				
Visual Question Answering	Visual Question Answering	Responds to various questions about the content of a radiology image.				
	Difference Visual Question Answering	Compares past (reference) and current (main) images to answer questions about the differences between them.				
	Visual Instruction-following	Generates responses by following specific instructions, facilitating clear communication in multiple interactions.				

Table 8: Task Description.

Task	Conversation Template
Single-image	User: radiology image: <image/> Which of the following findings are present in the radiology images? Findings: {findings } Assistant: {findings } User: Based on the previous conversation provide a description of the findings in the radiology image. Assistant: {report }
Multi-image	User: radiology images: { images } Which of the following findings are present in the radiology images? Findings: { findings } Assistant: { findings } User: Based on the previous conversation provide a description of the findings in the current follow-up radiology images. Assistant: { report }
Multi-study	User: prior radiology images: { prior images } prior radiology report: { prior report } follow-up images: { follow-up images } The radiology studies are given in chronological order. Which of the following findings are present in the current follow-up radiology images? Findings: { findings } Assistant: { findings } User: Based on the previous conversation provide a description of the findings in the current follow-up radiology images. Assistant: { report }
Disease Classification	User: radiology image: <image/> Which of the following findings are present in the radiology image? Findings: {findings } Assistant: {findings }
Finding Grounding	User: radiology image: <image/> Is {finding} present in the radiology image? If so, provide the bounding box coordinates of the region. Assistant: {bbox }
Grounded Finding	User: radiology image: <image/> Provide a finding name for this region. {bbox} Assistant: {finding }
Abnormality Detection	User: radiology image: <image/> Provide the bounding box coordinates of abnormal regions in the radiology image. Assistant: { bbox }
Multi Finding Grounding	User: radiology image: <image/> Which of the following findings are present in the radiology image? Provide the bounding box coordinates if present. Findings: {findings} Assistant: {findings, bboxes }
Organ Grounding	User: radiology image: <image/> Provide the bounding box coordinates of {organ} in the radiology image. Assistant: {bbox}
Grounded Organ	User: radiology image: <image/> Provide an organ name for this region. {bbox} Assistant: {organ }
Grounded Phrase Generation	User: radiology image: <image/> Provide a radiology report phrase for the region. {bbox} Assistant: {phrase }
Phrase Grounding	User: radiology image: <image/> Provide the bounding box coordinate of the region this phrase describes: { phrase } Assistant: { bbox }
Anatomical Region Grounding	User: radiology image: <image/> Provide the bounding box coordinate of the anatomical region. { name } Assistant: { bbox }
Grounded Anatomical Region	User: radiology image: <image/> Provide an anatomical region name for this region. {bbox} Assistant: {name}
Visual Question Answering	User: radiology image: <image/> Answer the question. {question} Assistant: {answer}
Difference Visual Question Answering	User: reference: <image/> main: <image/> Using the provided reference and main radiology images answer the following question. { question } Assistant: { answer }
Visual Instruction-following	User: radiology image: <image/> {question} Assistant: {answer}

Table 9: Conversation template for each task. Phrases inside brackets $\{\}$ are replaced with appropriate text based on data instances, with red-colored text indicating target outputs. < **image**> indicates the positions where visual tokens are inserted.

Task type	Task	Dataset	Train	Dataset ratio	Task type ratio	
Medical	Single-image	MIMIC-CXR	270,236	200.0		
Report	Multi-image	MIMIC-CXR	151,606	120.0	0.54	
Generation	Multi-study	MIMIC-CXR	68,373	60.0	_	
Image Understanding	Disease Classification	BRAX CheXpert VinDr-CXR ChestX-ray14 ChestX-Det10 SIIM RSNA COVID-19 Radiography	40,965 223,414 15,000 86,523 3,578 10,675 14,863 15,153	7.0 30.0 3.0 10.0 1.0 4.0 4.0 3.0		
	Finding Grounding	VinDr-CXR ChestX-ray14 ChestX-Det10 JSRT SIIM RSNA COVID-QU-Ex QaTa-COV19	15,000 50,500 3,578 247 10,675 14,863 2,796 4,194	2.0 10.0 2.0 0.1 4.0 1.0 1.0 1.0	_	
	Grounded Finding	VinDr-CXR ChestX-Det10 JSRT SIIM RSNA COVID-QU-Ex QaTa-COV19	4,394 2,967 154 2,379 6,012 1,864 4,194	2.0 2.0 0.1 2.0 1.0 1.0 3.0	0.35	
	Abnormality Detection	VinDr-CXR ChestX-Det10 JSRT SIIM RSNA COVID-QU-Ex QaTa-COV19	4,394 2,967 154 2,379 6,012 1,864 4,194	3.0 1.5 0.1 2.0 1.0 1.5 1.5		
	Multi Finding Grounding	ChestX-ray14 VinDr-CXR ChestX-Det10	50,500 15,000 3,578	10.0 2.0 2.0		
	Organ Grounding	COVID-19 Radiography COVID-QU-Ex	15,153 3,728	0.8 0.8		
	Grounded Organ	COVID-19 Radiography COVID-QU-Ex	15,153 3,728	0.8		
	Grounded Phrase Generation	MS-CXR ImaGanoma	638	2.0		
	Generation	MS-CXR	638	2.0	-	
	Phrase Grounding	ImaGenome	164,229	40.0		
	Anatomical Region Grounding	ImaGenome	164,229	20.0		
	Grounded Anatomical Region	ImaGenome	164,229	20.0		
Visual	Visual Question Answering	MIMIC-CXR-VQA MIMIC-Diff-VQA	255,919 553,156	2.0 2.0		
Question	Difference Visual Question Answering	MIMIC-Diff-VQA	129,900	4.0	0.11	
Answering	Visual Instruction-following	RaDialog	297,964	70.0	-	

Table 10: The number of training instances and sampling ratios for each dataset used across tasks.

Ground Truth

CheXbert positive labels: Atelectasis, pleural effusion, support devices.

FINDINGS: A left central venous catheter is seen terminating in the lower SVC. Again seen is elevation of the right hemidiaphragm and small pleural effusion and atelectasis at the base of the right lung. The cardiomediastinal silhouette and hilar contours are grossly unchanged. There is no evidence of pneumothorax. Thoracolumbar fusion hardware is seen unchanged in appearance.



Figure 5: Examples of medical report generation across various scenarios. For the same study, the top left shows the result for *single-image*, the top right for *multi-image*, and the bottom for *multi-study* report generation. The results are selected from the MIMIC-CXR test set.



Ground-truth Category: Atelectasis Text: there is atelectasis of the right upper lobe which may be due to bronchial obstruction and obstructing lesion **Bounding box**: [17,17,41,44]

Ground-truth



Pneumothorax.



Figure 6: Examples of visual grounding. The ground-truth bounding box is represented by a yellow solid box, while the prediction is shown with a red dashed box. The images are selected from the MS-CXR test set.



Ground-Truth

CheXbert positive labels: Lung Opacity FINDINGS: The heart is normal in size. The mediastinal and hilar contours appear within normal limits. Pleural effusions have resolved. There is a patchy new opacity in the lateral

segment of the right middle lobe, worrisome for pneumonia. However, elsewhere, the lungs appear clear. The osseous structures are unremarkable.



Figure 7: Examples of medical report explanations using easy language. The left shows the results from M4CXR, while the right shows the results when trained without VQA datasets.

M4CXR		Groun CheXb Suppor FINDI and is r pneumor remains	nd-Truth vert positive labels: Cardiomegaly, Pleural Effusion, t Devices. NGS: A right pleural effusion has increased since now large. The left lung is clear. No left effusion or othorax is present. A right-sided Port-A-Cath tip s in the mid SVC. Cardiomegaly is unchanged. v/o VQA
	Which of the following findings are present in the radiology image? Findings: Enlarged Cardiomediastinum,		Which of the following findings are present in the radiology image? Findings: Enlarged Cardiomediastinum,
	Cardiomegaly, Pleural Effusion, Support Devices.		Cardiomegaly, Pleural Effusion, Support Devices.
M4CAR	Based on the previous conversation, provide a description of the findings in the radiology image.	M4CAR	Based on the previous conversation, provide a description of the findings in the radiology image.
M4CXR	There is been substantial increase in the right pleural effusion since the prior radiograph. This however was already a substantial effusion in the prior study and may indicate congestive heart failure. Increased cardiac silhouette appears partially obscured by surrounding fluid. A right Port-A-Cath ends in the lower SVC. There is no left pleural effusion or left lung consolidation pulmonary edema or pneumothorax.	M4CXR	In comparison with the study of there is some increase in the degree of fluid within the right mid-to-lower zone with little change in the aeration of the right upper lobe. The central catheter again extends to the mid to lower portion of the SVC. Continued enlargement of the cardiac silhouette with little evidence of vascular congestion.
	Summarize the description in one concise sentence.		Summarize the description in one concise sentence.
()	There is a significant increase in the right pleural effusion, with cardiomegaly and the presence of a right Port-A-Cath.	()	Large right pleural effusion.
M4CXR		M4CXR	

Figure 8: Examples of medical report summarization. The left shows the results from M4CXR, while the right shows the results when trained without VQA datasets.



Ground-Truth

CheXbert positive labels: Atelectasis, Pleural Effusion. FINDINGS: In comparison with the study of ____, there is slightly less opacification at the left base. This is most consistent with atelectasis and effusion, though in the appropriate clinical setting, supervening pneumonia would have to be considered. Specifically, no evidence of pneumothorax. No acute focal pneumonia.



Figure 9: Examples of medical treatment recommendation. The left shows the results from M4CXR, while the right shows the results when trained without VQA datasets.